

# N3100 Series Nitrile Exam Gloves

## Features and Benefits

- Nitrile
- Thickness (Finger): 3.5mil
- Powder Free
- ASTM D6319, D6978
- Exam Grade
- Color: Blue
- Texture: Fingertip
- Length: Standard (240mm)
- Cuff: Beaded
- Non-Sterile
- Chemotherapy Tested<sup>1</sup>



<sup>1</sup>Chemotherapy drug testing information on back of box.

## Usage

Dependable for house-wide applications within Healthcare, Labs, EMS and other settings where high-use and disposable exam glove protection is important.

### Specifications



### Usage



First Responders



Healthcare



Industrial



Laboratory



Long Term Care

## Ordering Information

ITEM #	SIZE	GLOVES PER BOX	BOXES PER CASE	GLOVES PER CASE
N3100XS	X-Small	100	10	1000
N3100S	Small	100	10	1000
N3100M	Medium	100	10	1000
N3100L	Large	100	10	1000
N3100XL	X-Large	100	10	1000

# N3100 Series Nitrile Exam Gloves

100 CT



## Spec Information

OVERVIEW	
Material	Nitrile
Powder Content	Powder-Free
Freedom From Holes	1.5 AQL
Surface	Fingertip Textured
Weight (Based on Size M)	3.4 gram
Color	Blue
Cuff Length	240mm, Standard
Not Made With Natural Rubber Latex	Yes
Chemo Rated (Yes/No)	Yes, in accordance with FDA ASTM D6978

THICKNESS MEASUREMENTS	
Finger	.09mm / 3.5mm
Palm	.06mm / 2.5mm

PHYSICAL PROPERTIES		
Criteria	Before Aging	After Aging
Tensile Strength	>/= 14	>/= 14
Elongation	500 (min.)	400 (min.)

DIMENSIONS			
Item#	Size	Palm Width (mm)	Length
N3100XS	Extra-Small	70-80	240mm / 9.5 inches
N3100S	Small	81-90	240mm / 9.5 inches
N3100M	Medium	91-100	240mm / 9.5 inches
N3100L	Large	101-110	240mm / 9.5 inches
N3100XL	Extra-Large	111-120	240mm / 9.5 inches