



SH Medical Cabinets

School Health's modular medical cabinetry system is designed to complement our recovery couches and exam tables. Each cabinet features fully concealed, self-closing, three-way adjustable euro hinges and commercial-grade drawer glides. Security locks are standard for all doors and drawers.

Start off by choosing any combination of cabinetry from three expedient widths (19½", 25½", and 31½"), then simply choose one of our color-coordinated "E-Z Order Combos" (as illustrated below) or create your own combo. School Health now offers wire pull handles for all of our cabinetry and couches.

Refer to page 227 for available color and style options. To assure submission of all necessary info, use our Easy Cabinetry Design Order Form on page 269. The form will help you get exactly what you want, as quickly as possible.

E-Z Order Combos

Refer to the E-Z Combo number of your choice on your order. All E-Z combos are listed as follows:

- Cabinetry Color
- Countertop Color
- Pull Style



E-Z Combo 1

- Maple
- White
- Aluminum



E-Z Combo 2

- Custom Gray
- Fog Dust
- Aluminum



E-Z Combo 3

- Storm Gray
- Fog Dust
- Solid Maple



E-Z Combo 4

- Walnut
- White
- Wire Pulls



E-Z Combo 5

- White
- White
- White



E-Z Combo 6

- Almond
- Almond
- Solid Oak

* Direct-Ship Item. Please specify color when ordering; allow 3-6 weeks for fabrication and delivery after all information has been submitted. See index section for details.



25½" Cabinets

Colors shown are **E-Z Combo 2:**

- Cabinet Custom Gray
- Countertop Fog Dust
- Pulls Aluminum

E. 25½" Wall Cabinet *AC24098

N. 25½" Base Cabinet *AC24223



19½" Cabinets

Colors shown are **E-Z Combo 4:**

- Cabinet Walnut
- Countertop White
- Pulls Wire Pulls

A. 19½" Wall Cabinet *AC24209

M. 19½" Base Cabinet *AC24100



31½" Cabinets

Colors shown are **E-Z Combo 5:**

- Cabinet White
- Countertop White
- Pulls White

F. 31½" Wall Cabinet *AC24099

S. 31½" Sink Base with Towel Drop ... *AC24226



25½" Cabinets

Colors shown are **E-Z Combo 3:**

- Cabinet Storm Gray
- Countertop Fog Dust
- Pulls Solid Maple

K. 25½" Base Cabinet – 2 Drawers ... *AC24102



Wheelchair Accessible Sink Cabinet

Our wheelchair sink offers flexibility and can be mounted at any height, allowing for the recommended 27" of vertical knee clearance. Coordinates with any cabinet configuration. Standard removable protective panel underneath conceals plumbing fixtures.

Colors shown are **E-Z Combo I:**

Cabinet Maple
 Countertop White

T. 19 1/2" *AC24230



25 1/2" Cabinets

Colors shown:

Cabinet Custom Gray
 Countertop Custom Gray
 Pulls Wire Pulls

B. Wall Cabinet - Acrylic *AC24210

N. Base Cabinet - 2 Doors *AC24223

H. Base Cabinet - 5 Drawers *AC24221

SH Wall Cabinet Guide

**30"H x 12 1/4"D
19 1/2" Wide**



A. Clear Acrylic Framed Door
Wall Cabinet - Acrylic

*AC24209 \$398.00 ea.



D. Solid Laminate Door
Wall Cabinet - Solid

*AC24097 \$313.00 ea.

**30"H x 12 1/4"D
25 1/2" Wide**



B. Clear Acrylic Framed Doors
Wall Cabinet - Acrylic

*AC24210 \$499.00 ea.



E. Solid Laminate Doors
Wall Cabinet - Solid

*AC24098 \$382.00 ea.

**30"H x 12 1/4"D
31 1/2" Wide**



C. Clear Acrylic Framed Doors
Wall Cabinet - Acrylic

*AC24211 \$561.00 ea.



F. Solid Laminate Doors
Wall Cabinet - Solid

*AC24099 \$434.00 ea.

SH Base Cabinet Guide

**36"H x 19"D
19 1/2" Wide**




G. 5 Drawers

*AC24101 \$694.00 ea.



J. 1 Deep Drawer, 1 Door, 1 Adjustable Shelf

*AC24220 \$642.00 ea.



M. 1 Door, 2 Adjustable Shelves

*AC24100 \$520.00 ea.



T. Wheelchair Accessible Sink Cabinet

*AC24230 \$662.00 ea.

**36"H x 19"D
25 1/2" Wide**



H. 5 Drawers

*AC24221 \$784.00 ea.




K. 2 Doors, 2 Drawers, 1 Adjustable Shelf

*AC24102 \$749.00 ea.



N. 2 Doors, 2 Adjustable Shelves

*AC24223 \$599.00 ea.



P. 2 Doors, 1 Sink

*AC24227 \$883.00 ea.



R. 2 Doors, 1 Sink, 2 Towel Drop Doors

*AC24225 \$928.00 ea.

**36"H x 19"D
31 1/2" Wide**




I. 5 Drawers

*AC24222 \$882.00 ea.




L. 3 Drawers, 2 Doors

*AC24103 \$879.00 ea.



O. 2 Doors, 2 Adjustable Shelves

*AC24224 \$648.00 ea.



Q. 2 Doors, 1 Sink

*AC24228 \$917.00 ea.



S. 2 Doors, 1 Sink, 2 Dual Towel Drop Doors

*AC24226 \$969.00 ea.

These measurements are for standalone cabinets only.

* Direct-Ship Item. Please specify color when ordering; allow 3-6 weeks for fabrication and delivery after all information has been submitted. See index section for details.




NEW!
Exam Room Cabinets

Preconfigured cabinet groups allow reduction in cost by using only the components required and the options you choose. To offer better value these groups are offered in White, Almond, Folkstone, or Fog colors only. Panels, doors, and drawer fronts are edged in matching color 3mm PVC with rounded profile and corner rounding.

Options available are a stainless steel sink with goose neck faucet and strainer basket available on 24" and 30" cabinets only; work station slide out available on 18" and 24" cabinets only; and drawer or door locks (2 door cabinets require only 1 lock).





When ordering please submit item number, color and options, and position – left, right, or centered.


Cabinets come standard with wire pull handles.


A

B

C

F

G

H

D

E

K

I

J

L
Furnishings
Wall Cabinets – 30" H x 12" D

- A.** 18"W, 1 Door, and 2 Shelves **AC24881** \$245.00 ea.
- B.** 24"W, 2 Doors, and 2 Shelves **AC24882** \$288.00 ea.
- C.** 30"W, 2 Doors, and 2 Shelves **AC24883** \$325.00 ea.
- D.** 48"W, 3 Doors, and 4 Shelves **AC24884** \$479.00 ea.
- E.** 72"W, 5 Doors, and 6 Shelves **AC24885** \$730.00 ea.

Base Cabinets – 36" H x 18" D

- F.** 26"W, 1 Drawer, 2 Doors, and 1 Shelf **AC24886** \$545.00 ea.
- G.** 32"W, 1 Drawer, 2 Doors, and 1 Shelf **AC24887** \$645.00 ea.
- H.** 44"W, 2 Drawers, 3 Doors, and 2 Shelves **AC24888** \$929.00 ea.

50"W, 2 Drawers, 3 Doors, and 2 Shelves

- **AC24889** \$995.00 ea.
- I.** 8"W, 3 Drawers, 4 Doors, and 3 Shelves **AC24891** \$1,365.00 ea.
- J.** 74"W, 3 Drawers, 5 Doors, and 3 Shelves **AC24892** \$1,495.00 ea.

Available Options – Must be Ordered with Cabinets to Receive Price

- Lock** (2 Door cabinets Require Only 1 Lock) **AC24895** \$45.00 ea.
 - K.** Work Station Slide Out (For 18" and 24" Base Cabinets Only) **AC24894** \$69.00 ea.
 - L.** Stainless Steel Gooseneck Faucet and Strainer Basket (For 24" and 30" Base Cabinets only) ***AC24893** \$275.00 ea.
- * With the addition of a sink, one drawer will be removed from base cabinet.*

Exam Room Cabinet Colors




SH Corner Cabinets

Both the wall and base cabinets measure 36" along each back wall. One side of the cabinet is similar to cabinet M (see previous page) with two adjustable shelves, while the other side contains a handy wire "lazy susan" rack.




**30"H x 12 1/4"D
36" x 36" Wide**



U. Clear Acrylic Framed Doors

Corner Wall Cabinet

..... *AC24407 \$884.00 ea.



W. Solid Laminate Door

Corner Wall Cabinet

..... *AC24406 \$649.00 ea.

**36"H x 19"D
36" x 36" Wide**



**V. 2 Doors,
1 Lazy Susan,
2 Adjustable Shelves**

Corner Base Cabinet

..... *AC24408 \$915.00 ea.

SH Modular Cabinetry is custom-made and shipped directly to you from the manufacturer. Please allow 3-6 weeks (after all information has been submitted) for construction and delivery. To assure submission of all necessary information, please use our Cabinetry Design Order Form on page 269. The form will help you get exactly what you want, as quickly as possible!

E-Z Order Combos

Refer to the E-Z Combo number of your choice on your order. All E-Z combos are listed as follows:

- Cabinetry Color • Countertop Color • Pull Style



E-Z Combo 1

- Maple
- White
- Aluminum



E-Z Combo 2

- Custom Gray
- Fog Dust
- Aluminum



E-Z Combo 3

- Storm Gray
- Fog Dust
- Solid Maple



E-Z Combo 4

- Walnut
- White
- Wire Pulls



E-Z Combo 5

- White
- White
- White



E-Z Combo 6

- Almond
- Almond
- Solid Oak

Standard Custom


Cabinetry and Countertop Colors

These are the Standard Custom laminate colors School Health offers at no extra charge. Create your own Combo! Choose one each: Cabinetry, Countertop and Pull Style.

Cabinetry Colors



Countertop Colors



Pull Styles



LOOKING FOR FACILITY MANAGEMENT PRODUCTS?

Go to our **NEW Facility Management Section** starting on page 253!

* *Direct-Ship Item. Please specify color when ordering; allow 3-6 weeks for fabrication and delivery after all information has been submitted. See page 266 for details.*



SH Recovery Couches

Our Recovery Couches feature hardwood, dual-frame construction with separate base-headrest frames. Inner spring construction resists sagging, and is available on all but the Economy models. Paper roll dispensers, band cutters and exam paper included.

Recovery Couch Guide		Dimensions			Base Styles	
Type	Headrest Style	L	W	H	Solid Wood or Tubular Chrome legs	Deluxe Formica Base
Economy	Separate Pillow Included	72"	24"	19½"	✓ This Page	
Pediatric	Built-in	52"	27"	20"	✓ Page 232	
Varsity	Built-in	72"	27"	20"	✓ This Page	✓ Page 230
Biltmore	Adjustable	72"	27"	20"	✓ This Page	✓ Page 229
Race Car	Built-in or Adjustable	72"	27"	21"		✓ Page 232
Student	Separate Pillow Included	72"	29"	20"	✓ This Page	
Biltmore XL	Adjustable	76"	29"	20"	✓ This Page	
Palmer	Adjustable	76"	29"	23"	✓ This Page	




Varsity Express-Ship Colors

The Varsity Recovery Couch is our exclusive "Express Ship" Model. Select one of the "Express - Ship" colors, and your Varsity Couch will be shipped within 3-5 business days.



Palmer Recovery Couch

Featuring a larger headrest with dual-control locking device. This deluxe headrest includes a 12-position, fingertip control and positive locking device.
 Dimensions: 76"L x 29"W x 23"H
 Shipping Weight: 104 lbs.
 Palmer with Chrome Legs ***AC24137** \$483.00 ea.
 Palmer with Solid Wood Legs ***AC24046** \$495.00 ea.
 Clear Protective Couch Cover 76"L x 29"W **AC90578** \$28.25 ea.



Student Recovery Couch



Varsity Recovery Couch

Built-in headrest adds 6½" in height over body cushion.
 Dimensions: 72"L x 27"W x 20"H
 Shipping Weight: 95 lbs.

Varsity with Chrome Legs ***AC24048** \$435.00 ea.
 Varsity with Solid Wood Legs ***AC24142** \$448.00 ea.
 Clear Protective Couch Cover 72" L x 27" W **AC90504** \$28.25 ea.



B

Biltmore Recovery Couch

Both the Biltmore and Biltmore XL feature an adjustable 12-position headrest, fingertip control and positive locking device.

Biltmore Standard Recovery Couch

Dimensions: 72"L x 27"W x 20"H
 Shipping Weight: 95 lbs.
A. Standard with Chrome Legs ***AC24047** \$456.00 ea.
B. Standard with Solid Wood Legs ***AC24139** \$464.00 ea.
 Clear Protective Couch Cover 72"L x 27"W **AC90505** \$28.25 ea.

Biltmore XL Recovery Couch

Dimensions: 76"L x 29"W x 20"H
 Shipping Weight: 104 lbs.
 XL with Chrome Legs ***AC24250** \$484.00 ea.
 XL with Solid Wood Legs ***AC24251** \$495.00 ea.

Student Recovery Couch (Shown left)

One-inch square chrome-plated or solid wood legs securely mounted to a 4½" hardwood frame. Pillow included.
 Dimensions: 72"L x 29"W x 20"H
 Shipping Weight: 95 lbs.
 Student with Chrome Legs ***AC24051** \$348.00 ea.
 Student with Solid Wood Legs ***AC24138** \$372.00 ea.
 Clear Protective Couch Cover 72"L x 29"W **AC90519** \$28.25 ea.

Economy Couch (Not shown)

Same as Student Recovery Couch, with smaller dimensions.
 Dimensions: 72"L x 24"W x 19½"H
 Economy with Chrome Legs ***AC24052** \$319.00 ea.
 Economy with Solid Wood Legs ***AC24136** \$325.00 ea.

How To Order Recovery Couches

1. Select Leg Style



Solid Wood Legs



Chrome Legs

2. Select Upholstery

Naugahyde®

Decorative, vinyl coated fabric.



* Direct-Ship Item. Please specify color when ordering; allow 3-6 weeks for fabrication and delivery after all information has been submitted. See page 266 for details.

SIGNATURE RECOVERY COUCHES



A. Shown in Dutch Blue on Walnut Base with Wire Pulls

Built-In Concealed Compartment for Paper Roll Dispenser

B. Shown in Yew Green on Gray Base with Gray Pulls

C. Shown in Chianti Classic Twill (discontinued color) on Maple Base with Maple Pulls

D. Shown in Claret on Almond Base with Aluminum Pulls

Large, Lockable Drawers with Easy-Glide Rollers

SH Biltmore Custom Recovery Couch

Functionality meets design in our Custom Recovery Couches! Made with a durable, maintenance-free laminate base and outfitted with two, easy-glide drawers (20 3/8" W x 14 1/4" D x 4 1/2" H) and 12-position headrests.

Dimensions: 72"L x 27"W x 25"H

- A.** with Wire Drawer Pulls *AC24283 \$740.00 ea.
- B.** with Plastic Drawer Pulls (specify Gray, Almond or White) *AC24271 \$740.00 ea.

- C.** with Solid Wood Drawer Pulls (specify Oak or Maple) *AC24272 \$740.00 ea.
- D.** with Aluminum Drawer Pulls without Drawers (Not shown) *AC24270 \$740.00 ea.
- *AC24273 \$599.00 ea.

* Direct-Ship Item. Please specify color when ordering; allow 3-6 weeks for fabrication and delivery after all information has been submitted. See page 266 for details.



SIGNATURE RECOVERY COUCHES

How To Order Recovery Couches

1. Select Drawer Pull Style

If you are selecting a couch with drawers, there are several drawer pull styles available:



Be sure to specify desired pull style, upholstery color and base color when placing your order. Please allow 3-6 weeks for fabrication and delivery after all necessary information has been submitted.

2. Select Upholstery Color

Naugahyde® Selection

Decorative, vinyl coated fabric.



3. Select Base Color



A. Varsity shown in Forest on Almond Base with Aluminum Pulls



D. Varsity shown in Black on Almond Base without drawers



B. Varsity shown in Tea Rose on White Base with White Plastic Pulls



C. Varsity shown in Regimental Blue on Custom Gray Base with Oak Pulls

Furnishings

SH Varsity Custom Recovery Couch

School Health Varsity Recovery Couches are made with laminate bases and have two easy-glide drawers (locks included).

Dimensions: 72"L x 27"W x 25"H

- A.** with Aluminum Drawer Pulls ***AC24274** \$699.00 ea.
- B.** with Plastic Drawer Pulls (specify Gray, Almond or White) ***AC24275** \$699.00 ea.
- C.** with Solid Wood Drawer Pulls (specify Oak or Maple) ***AC24276** \$699.00 ea.
- with Wire Drawer Pulls (Not shown) ***AC24282** \$699.00 ea.
- D.** without Drawers ***AC24277** \$540.00 ea.

* Direct-Ship Item. Please specify color when ordering; allow 3-6 weeks for fabrication and delivery after all information has been submitted. See page 266 for details.

Recovery couches and treatment tables that fit your needs and your budget!


New and Improved!
Couches & Tables now come with a paper roll dispenser and have an upgraded, anti-microbial vinyl that is stain, mildew, and oil resistant.



Removable head rest fits securely with Velcro® straps.



A Shown with Navy (24817) upholstery option and wooden legs. Also available with chrome legs.



B Shown with Blue (24829) upholstery option and adjustable backrest. Convenient storage shelf.



C Shown with Gray (24833) upholstery option. Features adjustable back and legs.



Available Upholstery Colors



A. Recovery Couch

Features:

- Available in Chrome or Hardwood Legs
- 1" solid beech wood frame construction
- Woven, knit-backed vinyl upholstery
- 2" high-density foam to support the patient
- 350 lbs. weight capacity
- Shipping weight: 79 lbs.

72"L x 20"H x 26"W

Chrome Legs

Burgundy	*AC24820	\$315.00 ea.
Blue	*AC24824	\$315.00 ea.
Navy Blue.	*AC24822	\$315.00 ea.
Beige	*AC24823	\$315.00 ea.
Gray	*AC24821	\$315.00 ea.

Wooden Legs

Burgundy	*AC24815	\$315.00 ea.
Blue	*AC24819	\$315.00 ea.
Navy Blue.	*AC24817	\$315.00 ea.
Beige	*AC24818	\$315.00 ea.
Gray	*AC24816	\$315.00 ea.

B. Standard Treatment Table with Shelf

Features:

- Adjustable backrest – standard
- 2" high-density foam to support the patient
- Solid natural finish frame and legs
- Woven, knit-backed vinyl upholstery
- 350 lbs. weight capacity
- Shipping weight: 86 lbs.

74"L x 24"W x 28"H

Burgundy	*AC24825	\$445.00 ea.
Black	*AC24826	\$445.00 ea.
Blue	*AC24829	\$445.00 ea.
Gray	*AC24828	\$445.00 ea.
Navy Blue.	*AC24827	\$445.00 ea.

C. Adjustable Height Treatment Table

Features:

- Adjustable backrest – standard
- 2" high density foam to support the patient
- Solid natural finish frame and legs
- Woven, knit-backed vinyl upholstery
- 350 lbs. weight capacity
- Shipping weight: 86 lbs

74"L x 24"W x 23 1/2" – 35"H

Burgundy	*AC24830	\$445.00 ea.
Black	*AC24831	\$445.00 ea.
Blue	*AC24834	\$445.00 ea.
Gray	*AC24833	\$445.00 ea.
Navy Blue.	*AC24832	\$445.00 ea.

* Available for shipment on 4/30/2011 and will be delivered within 4-7 business days.



Furnishings



PEDIATRIC RECOVERY COUCHES & EXAM TABLES



SH Pediatric Recovery Couch

A smaller version of our popular Varsity model. Specify your color choice for upholstery when ordering.

Dimensions: 52"L x 27"W x 20"H

Shipping Weight: 75 lbs.

with Chrome Legs . . . *AC24252 \$339.00 ea.

with Hardwood Legs *AC24253 \$375.00 ea.

Clear Protective Couch Cover

52"L x 27"W AC90413 \$28.25 ea.



AC24082
Shown in Tomato



AC24039
Shown in Dutch Blue

Race Car Recovery Couch

Please specify color choice when ordering.

Dimensions: 72"L x 27"W x 21"H

Shipping Weight: 220 lbs.

with Built-In Headrest

..... *AC24082 \$689.00 ea.

with Adjustable 12-Position Headrest

..... *AC24039 \$699.00 ea.

Upholstery Color Choices



SH Pediatric Recovery Base Couch with Drawer

A smaller version of our popular base model couches. Please see page 230 for upholstery color and base choices.

Dimensions: 52"L x 27"W x 24"H

with Aluminum Pulls. *AC24295 \$669.00 ea.

with Solid Wood Pulls (specify Oak or Maple)

..... *AC24297 \$594.00 ea.

Solid Base without Drawer

..... *AC24299 \$565.00 ea.

Couch Extension Piece

Use for additional seating to coordinate with your health office or as an additional piece to add to your pediatric couch. See page 173 for available upholstery colors and base choices.

20"L x 27"W x 24"H

..... *AC24373 \$253.00 ea.



Fire Engine & School Bus Pediatric Exam Tables

Fire Engine Treatment Table

Dimensions: 68"L x 24"W x 34"H

Shipping Weight: 350 lbs.

with 5 Drawers and Cabinet Storage

..... *AC24040 \$1,391.00 ea.

with 5 Drawers, Cabinet Storage, and Optional

Backrest (Not shown)

..... *AC24090 \$1,524.00 ea.

School Bus Treatment Table

Dimensions: 68"L x 24"W x 34"H

Shipping Weight: 350 lbs.

with 5 Drawers and Cabinet Storage

..... *AC24050 \$1,185.00 ea.

* Direct-Ship Item. Please specify color when ordering; allow 3-6 weeks for fabrication and delivery after all information has been submitted. See page 266 for details.

NEW!

Introducing Green-Line Eco-Friendly Products

An environmentally safer alternative to traditional materials for improved air quality.

- Exceeds CARB I and II standards
- All Green-Line wood has no added urea-formaldehyde
- All urethane fabric is CAL 01350
- PVC-Free construction.

Eco-Friendly Treatment Tables

- H-Brace is recessed into table leg with unique 4-sided "Lock-Tite" joint.
- Fusion Maple high pressure laminate legs, apron and stretchers.

Treatment Table

- 72" L X 24" W X 31" H
 *AC24835 \$425.00 ea.
 72" L X 30" W X 31" H
 *AC24836 \$470.00 ea.

Treatment Table with Shelf

- 72" L X 24" W X 31" H
 *AC24837 \$499.00 ea.
 72" L X 30" W X 31" H
 *AC24838 \$549.00 ea.

Treatment Table with Cabinet

- Fully-enclosed, all-laminate cabinet with sliding laminate doors.
 - Interior has divider and (1) adjustable shelf.
 - Fusion Maple Laminate.
- 72" L X 30" W X 31" H
 *AC24839 \$1,130.00 ea.

Accessories

Optional backrest and paper dispenser are not included. Backrest must be ordered at the same time as the treatment tables.

- A. *Optional backrest w/ gas springs*
 AC24896 \$185.50 ea.
 B. *Optional paper dispenser*
 AC24897 \$35.00 ea.

Eco-Friendly Carts

- (2) 5 1/2" Drawers & (1) 8 1/2" deep Drawer.
 - Maintenance-free Fusion Maple laminated surfaces on all cabinet exteriors and interiors.
 - 2 1/2" twin wheel hooded swivel Nylon Casters.
- 18.5" L X 18.5" W X 30" H
 *AC24841 \$475.00 ea.

Eco-Friendly Recovery Storage Couch

- Cabinet with (2) sliding laminate doors and finished back. Provides convenient storage space for all blankets, gowns, and supplies.
 - Folkstone Gray laminate finish.
- 72" L X 24" W X 22" H
 *AC24842 \$570.00 ea.

Accessories

Optional backrest and paper dispenser are not included. Headrest must be ordered at the same time as the recovery couch.

- C. *Optional headrest.* AC24898 \$149.75 ea.
 B. *Optional paper dispenser*
 AC24897 \$35.00 ea.



Shown with optional headrest.

PVC-FREE Urethane Top Color Options



Indigo Blue Flint Grey Mushroom Port Beige

* Direct-Ship Item. Please specify color when ordering; allow 3-6 weeks for fabrication and delivery after all information has been submitted. See page 266 for details.



PEDIATRIC EXAM TABLES



Toddler Changing Table and Wall Mounted Storage

Changing Table

The changing table features pull-out stairs for toddlers, which conveniently glides under table when not in use. Extra deep rails exceed recommended safety standards. Three generous supply shelves, one fixed and two adjustable. Paper roll dispenser included accommodates a fifteen-inch paper roll and attaches to either end of the table. Changing pad and safety strap included.

Dimensions: 22½" W x 48" L x 39" H; Weight: 145 lbs.
Changing Table *AC24543 \$780.50 ea.

Wall Storage

Install the convenient wall mounted storage for diapering necessities above changing table for a complete changing station. Mounting hardware included.

Dimensions: 48"W x 12"D x 25"H; Weight 64 lbs.
Wall Mounted Storage *AC24544 \$238.00 ea.

Combo Table and Storage

. *AC24545 \$988.00 ea.



Portable Sink Hand Washing Station

Use to encourage hand washing habits without the plumbing. Water is heated 20° F by the tankless, on-demand water heater. Faucet is ADA compliant. Sinks are 5" deep, made of FDA approved plastic. Other safety features include UL and NSF approved components, power switch inside the cabinet, coin-locked doors, four heavy duty locking casters and openings for cord management. One-year warranty. Choose from Birch wood or MapleWeave™ Laminate.

Birch, Child Sink

Size: 28" W x 23½" D x 24" H; Weight: 105 lbs.
. *AC24539 \$1,708.00 ea.

Birch, Adult Sink

Size: 28" W x 23½" D x 38" H; Weight: 135 lbs.
. *AC24540 \$1,743.00 ea.

MapleWeave Laminate, Child Sink

Size: 28" W x 23½" D x 24" H; Weight: 109 lbs.
. *AC24541 \$1,691.00 ea.

MapleWeave Laminate, Adult Sink

Size: 28" W x 23½" D x 38" H; Weight: 149 lbs.
. *AC24542 \$1,726.00 ea.



Pediatric Table

Solid wood construction. Table reaches a convenient treatment height of 36". Specify Natural or Walnut veneer finish as base color, and upholstery color.

Dimensions: 48"L x 24"W x 36"H

Shipping Weight: 120 lbs.
with Full-Length Shelf and 2 Fixed Storage Shelves (Shown)

. *AC24189 \$535.00 ea.
with Full-Length Shelf

. *AC24193 \$395.00 ea.
with Enclosed Central Cabinet and 2 Shelves

. *AC24190 \$575.00 ea.



Deluxe Pediatric Exam Table with Tanita Scale

Table with three-sided protective rim and 36" long measuring ruler. Scale has capacity of 44 lbs. Measures in 1 oz. increments from 0-20 lbs. and 2 oz. increments from 21-44 lbs. Performs self-check.

Dimensions: 48"L x 24" W x 36"H

Shipping Weight: 220 lbs.
. *AC24328 \$1,133.00 ea.



Deluxe Pediatric Exam Table

With one fixed shelf; no scale.

Dimensions: 48"L x 24" W x 36"H

Shipping Weight: 220 lbs.
. *AC24326 \$659.00 ea.

* Direct-Ship Item. Please specify color when ordering; allow 3-6 weeks for fabrication and delivery after all information has been submitted. See page 266 for details.

SH Pediatric Tables

Our pediatric exam tables offer the same quality construction as our custom models. They feature a high-density urethane foam top for maximum comfort and support, and the convenience of a built-in paper dispenser and cutter. Your choices of laminate and Naugahyde® upholstery are shown below.

How To Order Exam Tables

1. Select Drawer Pull Style



2. Select Upholstery

Naugahyde® Selection. Decorative vinyl coated fabric.



3. Select Base Color

Any above upholstery color can be ordered with any of our six available base colors (Maple, Oak, Walnut, Gray, Almond, White). See available custom base colors on page 230. Simply indicate your color choices on your order.



A. Rubbermaid® Sturdy Station 2 Changing Table

Can be mounted on a wall. Large, deep bed with adjustable safety belt for easy use and caution. Rated for use up to 50 lbs. Not a custom piece.

Size open: 35½"L x 19½"D x 28¼"W

Size closed: 35½"L x 4"D x 28¼"W

. AC24199 \$275.00 ea.

B. Protective Liners for Rubbermaid® Sturdy Station 2

Two-Ply Tissue Paper
320/case. AC21211 \$109.00 cs.

EXAM TABLES

Exam Tables

Our exam tables provide an array of designs and arrangements to suit your health room needs. They feature security locks, three-sided protective rims, spacious cabinets, and your choice of door and drawer pulls. The exam tables shown are finished in popular laminate colors. Locking systems are standard in all our cabinetry. Color choices are shown to the left.



Deluxe Space-Saver Exam Table

Specify pulls, upholstery, and base choice from chart on previous page.

Dimensions when open: 68"L x 24"W x 32"H
Length when closed: 51"

with Adjustable Footrest

.....*AC24181 \$629.00 ea.



Space-Saver Exam Table

Our standard Space-Saver Exam Table with fold-down back and footrests maximizes the smallest of areas. Features include a 2" foam top, full-length shelf, and upholstery and base in your choice (see page 236).

Dimensions when open: 68"L x 24"W x 30"H

Length when closed: 51"

with Adjustable Backrest and Footrest

.....*AC24179 \$556.00 ea.

with Adjustable Footrest

.....*AC24177 \$525.00 ea.



Adapta 300 Treatment Table

Features three highly flexible sections. Head section that moves from +45° to -90°, and mid-section ranges from 10° to 50° in 10° increments. Foot section ranges from 0° to 75° in six repeatable steps. Optional factory-installed casters available.

Three Sections: 28" x 75"

Adjustable Height: 19" to 37"

Shipping Weight: 250 lbs.

Midnight Grey*AC24344 \$2,595.00 ea.



Ritter® 204 Manual Exam Table

Offers efficient patient care. 18-gauge steel shell reinforced seat to support 500 lb. weight capacity and prevent breakage. Multi-density foam top design and comfort stirrups help make the lithotomy position a little easier for your patients. Powder-coated paint finish.

Dusty Blue, Moss, Navy, Perfect Plum, Shadow, Clay, Pebble Gray

.....*AC24360 \$1,335.00 ea.



Hi-Lo Changing Treatment Table

Ideal for accommodating children with special needs. One section top piece and side rails that pivot down for easy patient transfer and access; sturdy, powder-coated X-Frame design; moves smoothly up and down with no horizontal travel. Wheelchair accessible height range is 18" - 37". Gray vinyl with 1½" firm density foam and 4 swivel locking casters. "Hands-Free" Footswitch control operation.

Dimensions: 400 lb. capacity; 120-volt 60hz motor; 60"L or 72"L x 25"W x 18"H - 37"H

60" Length*AC24557 \$2,235.00 ea.

72" Length*AC24558 \$2,395.00 ea.

**Optional Belt-Strap

.....*AC24559 \$115.00 ea.

**Optional Belt-Strap must be ordered at time of purchase of table. This can not be added at a later time.


Furniture Accessories



Flat Pillow

Specify a color from the chart on page 236.

14" x 14"*AC24054 \$24.00 ea.



Wedge Pillow

Specify a color from the chart on page 236.

25" x 25" x 7", tapered to 1"
.....*AC24106 \$54.00 ea.

Additional Paper Roll Dispenser

with Plastic Band Cutter (for mounting under head)
.....*AC24049 \$24.00 ea.

Plastic Band Cutter Only
.....*AC24095 \$5.00 ea.

Clear Protective Couch Cover

Keep your couch looking like new!
12 gauge thickness
72"L x 24" W*AC90269 \$28.00 ea.

Fantastik® All Purpose Cleaner

Perfect for Twill, Naugahyde® and laminates.
32 oz.*AC34195 \$5.20 ea.

Vinyl Liquid Patch

1 oz.*AC90068 \$4.30 ea.

* Direct-Ship Item. Please specify color when ordering; allow 3-6 weeks for fabrication and delivery after all information has been submitted. See page 266 for details.



Features:

- Paper dispenser and paper cutter.
- Latex-free.
- 2" thick urethane foam top allows for maximum patient support and comfort. Foam conforms to patient then rebounds to original shape.
- All hardwood leg construction with wood veneer.

How To Order Treatment Tables

I. Select Base Veneer Finish



Natural



Walnut

Please Note: Colors are applied to a wood veneer and will differ slightly in appearance from a laminate of the same color.

2. Select Upholstery

Naughahyde®

Decorative, vinyl coated fabric.



Tea Rose



Mandarin Orange



Sun Yellow



Tomato



Graystone



Dove




Crocus



Claret



Adobe White



Black



Yew Green



Cinnabar



Sand



Royal



Bronze



Burgundy



Blue Ridge



Forest



Dutch Blue



Regimental Blue



Turquoise



Dusty Jade



Teal


* Direct-Ship Item. Please specify colors when ordering; allow 3-6 weeks for fabrication and delivery after all necessary information has been submitted. See page 266 for details.



SH Professional Treatment Table with Open Storage

All tables are 78" long, 30" high and have a full-length shelf.

Appearance	Available Options			SH Item #	Price
	Backrest	Face Cutout	Two Fixed Shelves		
A. with 24" Width	✓	✓	✓	*AC24212	\$514.00 ea.
with 30" Width	✓	✓	✓	*AC24191	\$589.00 ea.
B. with 24" Width		✓	✓	*AC24213	\$645.00 ea.
with 30" Width		✓	✓	*AC24192	\$719.00 ea.
C. with 30" Width	✓	✓		*AC24185	\$635.00 ea.
with 30" Width and Backrest		✓		*AC24186	\$775.00 ea.
D. with 24" Width	✓		✓	*AC24320	\$595.00 ea.
with 30" Width	✓		✓	*AC24321	\$657.00 ea.
with 24" Width and Backrest			✓	*AC24322	\$697.00 ea.
with 30" Width and Backrest			✓	*AC24324	\$778.00 ea.



SH Deluxe Professional Treatment Table with Enclosed Storage

Description	Length	Width	Height	SH Item #	Price
E. with Enclosed Cabinet, Sliding Doors and One Fixed-Position Shelf					
(Shown)	78"	30"	30"	*AC24315	\$1,135.00 ea.
Same as AC24315 with Backrest	78"	30"	30"	*AC24316	\$1,230.00 ea.
F. with Enclosed Central Cabinet and Two Shelves					
(Shown)	78"	30"	30"	*AC24187	\$840.00 ea.
Same as AC24187 with Backrest	78"	30"	30"	*AC24188	\$878.00 ea.
G. with Three 20" Wide Drawers on Nylon Rollers with Positive Stops					
(Shown)	78"	30"	30"	*AC24317	\$950.00 ea.
Same as AC24317 with Backrest	78"	30"	30"	*AC24318	\$1,025.00 ea.

Furnishings

NEW!

Electric Hospital Beds



- Features many positions such as Trendelenburg and Reverse Trendelenburg – ideal for many healthcare facilities.
- Electrostatic powder-coated antimicrobial (high concentrations of silver ions) bed surface provides longer lasting, durable finish
- Hospital Grade Grounded Plug is Standard Mattress Deck Length – 76"/80" optional 4" extender available
- Mattress Deck Width – 35"/optional 42" extender
- Overall Bed Length – 85.5" with bumper
- Overall Bed Width – 39.75" with rails/assist
- Maximum back angle – 70°
- Maximum knee angle – 30°
- Maximum multiple height – 30"
- Minimum height – 8.9"
- Horizontal translation – less than 1"
- Safe Workload (including resident, support surface, head/foot boards, and accessories) – 500 lbs.

- Strong Supportive grid deck is superior in frame design and offers a "no-sag" mattress support surface
- Patent pending design keeps bed level and balanced at all times providing a more comfortable, restful sleep
- 8.9 inch to 30 inch travel speed allows for bed positioning in as little as 22 seconds
- Straight vertical travel takes up less room space
- 360° mobility at ANY height – can move from room to room at any height even the lowest
- Adjustable leg lift system comes standard with beds for improving comfort for residents with edema and poor leg circulation
- Wall Saver is adjustable for close positioning to walls and keeps your walls scratch free
- Structured Base provides a reinforced, robust, base that is designed for support yet sufficient underbed clearance, ideal for using overbed tables and patient lifts
- Designed to support the use as a transport device – ideal for Sub-Acute and Rehab facility needs

4" Standard Casters (2 Locking) and Staff Control

76" Grid Deck	AC24855	\$1,643.00 ea.
80" Grid Deck	AC24856	\$1,643.00 ea.
76" Pan Deck.....	AC24857	\$1,727.00 ea.
80" Pan Deck.....	AC24858	\$1,727.00 ea.

One-Step Locking Pedals and Staff Control

76" Grid Deck	AC24589	\$1,799.00 ea.
80" Grid Deck	AC24561	\$1,799.00 ea.
76" Pan Deck.....	AC24562	\$1,883.00 ea.
80" Pan Deck.....	AC24563	\$1,883.00 ea.

4" Standard Casters (4 Locking) and Staff Control

76" Grid Deck	AC24864	\$1,664.00 ea.
80" Grid Deck	AC24865	\$1,664.00 ea.
76" Pan Deck.....	AC24866	\$1,750.00 ea.
80" Pan Deck.....	AC24867	\$1,750.00 ea.

One-Step Locking Pedals, Back-lit Hand Set, Under-bed Light and Staff Control

76" Grid Deck	AC24868	\$1,871.00 ea.
80" Grid Deck	AC24869	\$1,871.00 ea.
76" Pan Deck.....	AC24871	\$1,957.00 ea.
80" Pan Deck.....	AC24862	\$1,957.00 ea.

Foam Mattresses

76" L x 35" W.....	AC24873	\$255.00 ea.
80" L x 35" W.....	AC24874	\$255.00 ea.
76" L x 35" W with Zipper.....	AC24875	\$269.00 ea.
80" L x 35" W with Zipper.....	AC24876	\$269.00 ea.

Head and Foot Board Sets with Embedded Staff Control

Solar Oak, Rectangular Shaped	AC24877	\$249.00 set
Walnut, Rectangular Shaped	AC24878	\$249.00 set
Mahogany, Center Crown Shaped	AC24879	\$249.00 set
Cherry, Center Crown Shaped	AC24880	\$249.00 set

Electric Hospital Beds are custom made for each order, therefore School Health cannot accept returns.