




AAB GROUP

**AAB**

# DISPOSABLE PROTECTIVE MASK FOR CHILDREN

AAB (China) CO.,Ltd.

AAB 20 Years Experience





# ENTERPRISE QUALIFICATION

企业资质

# About AAB

'A-grade Exporting Enterprise'



✓ 'No.1 Hygiene Products Exporter' in China

AAB (China) CO. Ltd. is a professional manufacturer of hygiene products, including Baby diapers, baby pull-up, adult-diapers, incontinent pads, sanitary napkins, panty liners, tissue products, wet wipes, mask, dry and wet towels, with an area more than 80,000m<sup>2</sup>, and staffs more than 800.

AAB set up Global Joint R&D Center and spent 3 years to do resource integration from 8 departments of company & raw material suppliers for the purpose of offering one stop service for customer.

Our products are shipped to regions over Europe, America, South America, the Middle East, Southeast Asia and Africa with our OEM & ODM service and overall solutions.



# Company Registration



**营业执照**  
(副本) 副本编号: 1-1

  
扫描二维码登录  
“国家企业信用信  
息公示系统”了解  
更多登记、备案、  
许可、监管信息。

统一社会信用代码	91350583777549414P		
名称	中天(中国)工业有限公司	注册资本	23600.000000万港币
类型	有限责任公司(港澳台法人独资)	成立日期	2005年12月02日
法定代表人	黄家齐	营业期限	2005年12月02日至2055年12月01日
经营范围	生产生活用品兼分装,研发、生产、销售;医用防护口罩、非医用日用防护口罩;货物或技术的进出口业务(依法须经批准的项目,经相关部门批准后方可开展经营活动)	住所	南安市洪塘镇中天工业园
		登记机关	 2020年 3月 9日

国家企业信用信息公示系统网址: <http://www.gsxt.gov.cn>

市场主体应当于每年1月1日至6月30日通过  
国家企业信用信息公示系统报送公示年度报告

国家市场监督管理总局监制



# International Certificate

## FDA




2020

### CERTIFICATE OF REGISTRATION

This Certifies that

**AAB (CHINA) CO., LTD.**  
 AAB Industrial Park, Honglai Town, Nanan City, Fujian, 362331, China

Was registered with US Food & Drug Administration, Center for Devices and Radiological Health, pursuant to the Code of Federal Regulations 21 CFR 807. Such registration has been verified, with the registrant's authorization, by Manton Business and Technology Services.

**Registration Number:** 3008111543

Classification Name	Product Code	Product Name
Respirator, Surgical	MSH	6895 Respirator
Accessory, Surgical/Apparel	LYU	face mask
Mask, Scavenging	KHA	Isolation Mask
Pad, Menstrual, Unscented	HHDD	Sanitary Napkin
Pads, Menstrual, Scented-Deodorized	NRC	Panty Liner
Bedding, Disposable, Medical	KME	Under Pad
Fiber, Medical, Absorbent	FRL	Incontinent Pad
Mask, Oxygen	BYG	Oxygen Mask

**Date of Initial verification:** January 03, 2020  
**Date of Update:** April 01, 2020

*This Certificate affirms that Manton Business and Technology Services has verified that the above stated facility is registered with the US Food & Drug Administration, Center for Devices and Radiological Health, pursuant to the Code of Federal Regulations 21 CFR 807 on the date stated above, and makes no other representations and warranties, nor does this certificate makes other representations and warranties to other person or entity other than the name certificate holder, for whose sole benefit it is issued. Manton Business and Technology Services assume no liability to any person or entity in connection with the foregoing. Manton Business and Technology Services is a private registration agent and is not affiliated with the US Food and Drug Administration.*

Manton Business and Technology Services  
 New Jersey, USA  
[www.manton-tech.com](http://www.manton-tech.com)

*Charles Shen*  
 Charles Shen  
 Director



# International Certificate

## CE for disposable Protective mask

شهادة - 증명서 - Certificat - 證明書 - Сертификат - Certificate

Form QAT\_10-M04, version 00, effective since March 6th, 2020

### Certificate of Compliance



No. 0P200315.AABUT65

Technical Construction File no. TPFJ20031122194

**Certificate's Holder:** AAB (CHINA) CO.,LTD.  
 AAB Industrial Park, Honglai Town, Nan'an City  
 Fujian Province, China

**Certification ECM Mark:**

**Product Model:** Disposable protective mask  
 FFP1, FFP2, FFP3

**Verification to:** Standard:  
 EN 149: 2001+A1:2009  
 related to CE Directive(s):  
 R 2016/425 (Personal Protective Equipment)

**Remark:** This document has been issued on a voluntary basis and upon request of the manufacturer. It is our opinion that the technical documentation received from the manufacturer is satisfactory for the requirements of the ECM Certification Mark. The conformity mark above can be affixed on the products accordingly to the ECM regulation about its release and its use.

Additional information and clarification about the Marking:



The manufacturer is responsible for the CE Marking process; this document has been issued on the basis of the regulation on ECM Voluntary Mark for the certification of products. RG01\_ECM rev.3 available at: [www.entecerma.it](http://www.entecerma.it)

**Issuance date:** 15 March 2020

**Expiry date:** 14 March 2025

Reviewer  
 Technical expert  
 Amanda Payne



Approver  
 ECM Service Director  
 Luca Bedonni



**Ente Certificazione Macchine Srl**

Via Ca' Bella, 243 - Loc. Castello di Serravalle - 40053 Valsamoggia (BO) - ITALY  
 ☎ +39 051 6705141 ☎ +39 051 6705156 ✉ info@entecerma.it 🌐 www.entecerma.it

# International Certificate

## CE for disposable medical mask & KN95 protective mask

Form QAT\_10-M04, version 00, effective since March 6th, 2020

### CE Documentation Review



No. 0P200315.AABUT66

**Holder:** AAB (CHINA) CO., LTD.  
AAB Industrial Park, Honglai Town, Nan'an City,  
Fujian Province, China

**Review goal:** Verification of the presence of the  
Technical File in regards of the Medical  
Devices Directive 93/42/EEC Annex VII

**Product:** Disposable medical mask (no sterile)  
**Model(s):** KN95 mask, Medical disposable mask,  
Medical disposable surgical mask

**Classification:** Class I (no sterile)  
(accordingly to the Manufacturer's declaration)

**Review output:** We attest that a Technical File in reference to the  
Directive 93/42/EEC is in place for the CE Marking  
process. Technical File identified with the no.  
**TMFJ20031122193.**  
The manufacturer is responsible for the CE Marking  
process, and not exempted to carry out all necessary  
compliance activities. This document has been issued  
on the basis of the regulation on ECM Voluntary Mark  
for the certification of products. RG01\_ECM rev.3  
available at: [www.entecerma.it](http://www.entecerma.it)

**Issuance date:** 15 March 2020

**Expiry date:** 14 March 2025

Reviewer  
Technical expert  
Amanda Payne

Approver  
ECM Service Director  
Luca Bedonni

Review Report - 审查报告 - 검토 보고서 - Rapport d'Evaluation



# International Certificate

## SA8000

**BUREAU VERITAS**  
Certification



**AAB (CHINA) CO.,LTD.**

AAB Industrial Park, Honglai Town, Nan'an City, Fujian Province, China  
Bureau Veritas Certification certify that the Management System of the above organisation has been audited and found to be in accordance with the requirements of the management system standard detailed below  
*Standard*

**SA 8000:2014**

*Scope of certification*

**MANUFACTURE OF SANITARY NAPKINS, SANITARY PADS, DIAPER, PAPER PADS**

Certification cycle start date: **11 Sep 2018**

Subject to the continued satisfactory operation of the organisation's Management System, this certificate expires on: **10 Sep 2021**

Original certification date: **11 Sep 2018**

Certificate No. **IND18.10046** Version : 1 Revision date: **11 Sep 2018**

*Certification Authority*



*Local office Address*

F22, Tower B, Beijing Global Trade Center, 36 North Third Ring Road East, Dongcheng District, Beijing China, 100013

*Managing office address*

Bureau Veritas (India) Private Limited (Certification Business)  
72 Business Park, Marol Industrial Area, MIDC Cross Road 'C', Andheri (East), Mumbai - 400 059, India.

*Disclaimer :*

"Social Accountability International and other stakeholders in the SA8000 process only recognize SA8000 certificates issued by qualified CBs granted accreditation by SAAS and do not recognize the validity of SA8000 certificates issued by unaccredited organizations or organizations accredited by any entity other than SAAS."

To check this certificate validity please call: **+861059683663** or visit SAAS website @ [www.saasaccreditation.org/certification](http://www.saasaccreditation.org/certification)

Further clarifications regarding the scope of this certificate and the applicability of the management system requirements may be obtained by consulting the organisation.





# International Certificate

## GMP



# CERTIFICATE OF REGISTRATION

This is to certify that the quality management system of:

## AAB (CHINA) CO., LTD.

2/F and 3/F, Area B of Production Building, AAB Industrial Park, Honglai Town, Nan'an City, Fujian Province, China

has been assessed by Intertek as conforming to the requirements of:

**CURRENT GOOD MANUFACTURING PRACTICE BASED ON  
21 CFR PART 820 (2017) PUBLISHED BY U.S. FOOD AND  
DRUG ADMINISTRATION**

The scope of activities:

Manufacturing of sanitary napkin, panty liner and disposable diaper



**Intertek**

Certificate Number:

SZ1809E3

Initial Audit Date:

23-25 Jul, 2018

Certificate Issue Date:

28 Sep, 2018

Certificate Expiry Date:

27 Sep, 2021

Certification Administration Centre Intertek  
Testing Services



This certificate is valid as long as it bears a proper and authentic Intertek's Laser Logo dedicated for the year of initial certification and after satisfactory annual surveillance.

Helen Xue  
General Manager  
Chemicals & Pharmaceuticals Division





# PRODUCT DISPLAY

产品展示



# DISPOSABLE PROTECTIVE MASK FOR CHILDREN

AAB儿童平面口罩



Comfortable  
ear with



Non-woven fabric  
air permeability



3 ply  
protection



# DISPOSABLE PROTECTIVE MASK FOR CHILDREN


AAB儿童平面口罩





# ERGONOMIC TAILORING, COMFORTABLE WEARING EXPERIENCE

符合人体工学的裁剪工艺，佩戴舒适







## SMS STRUCTURE

### NON-WOVEN FABRIC LAYER


SMS结构 无纺布层





# MELT-BLOWN NON WOVEN CLOTH

熔喷布





## TWO COLOURS RANDOM DELIVERY

两种颜色随机发货





## OF THE GOODS

DISPOSABLE PROTECTIVE MASK

## PERIOD OF VALIDITY

TWO YEARS

## PRODUCT SIZE

XCM( $\pm 0.5$ CM)

## PACKING QTY.

50 PCS/BOX,40 BOXES/CARTON

## QUALITY STANDARD

Q/AABQ 14-2020

## QTY./CARTON

PCS

## G.W.

KG

## CARTON SIZE

17\*10\*10CM

## VOLUME

0.084CBM

## 20 GP

700,000 PCS

## 40 HQ

1,700,000 PCS





# Production process

Merge three layers








# Production process

Sewing the nose clip, folding





# Production process


Tablet





# Production process

Cutting edge









# Production process

Hot press fixed hanger on the rope





# Workshop show







# PRODUCT QUALITY INSPECTION AND EXPORT QUALIFICATIONS

产品质量检测与出口资质

# Test Report



## TEST REPORT

(Electronic version)



VERIFICATION WEBSITE: [www.gtto.net.cn](http://www.gtto.net.cn)

VERIFICATION CODE: FIVG-3466-04



No:200041161

ISSUE DATE:2020-09-20

APPLICANT: AAB(CHINA)CO.,LTD  
 ADDRESS: 2F, Building C06, Huachuang Industrial Park, Shijie Town, Panyu District, Guangzhou

APPLICANT PROVIDED SAMPLE DESCRIPTION:  
 FORTY(40) PIECES OF DISPOSABLE PROTECTIVE MASK  
 BRAND: AAB  
 COLOUR: BLUE  
 MANUFACTURE'S NAME: AAB(CHINA)CO.,LTD

DATE RECEIVED/DATE TEST STARTED: 2020-09-16

CONCLUSION:

pH VALUE	M
FORMALDEHYDE CONTENT	M
BACTERIAL FILTRATION EFFICIENCY	M
PARTICLE FILTRATION EFFICIENCY	M
INSPIRATORY RESISTANCE	M
EXPIRATORY RESISTANCE	M
APPEARANCE QUALITY[3 PIECES]	M

NOTE: "M"-MEET THE STANDARD'S REQUIREMENT "F"-FAIL TO MEET THE STANDARD'S REQUIREMENT  
 "----"-NO COMMENT

REMARK:

ALL THE TESTED ITEMS ARE TESTED UNDER THE STANDARD CONDITION (EXCEPT FOR INDICATION).  
 COPIES OF THE REPORT ARE VALID ONLY RE-STAMPED.  
 THE EXPERIMENT WAS CARRIED OUT AT No.1, CHUJIANG ROAD, PANYU DISTRICT, GUANGZHOU, GUANGDONG, P.R.CHINA.

APPROVED BY: 郭子山  
 Tashan Gao SENIOR ENGINEER

PAGE 1 OF 3

# Test Report



## TEST REPORT

(Electronic version)

No:200041161



PAGE 2 OF 3



## TEST REPORT

(Electronic version)

No:200041161

### pH VALUE

(GB/T 7873-2009, 0.1mol/L KCl EXTRACTING SOLUTION)  
(GB/T 7873-2009 MOD ISO 3071:2006)

6.3

REQUIREMENT  
4.0~8.5  
(Q/AADQ 15-2020)

### FORMALDEHYDE CONTENT

(GB/T 2912.1-2009, WATER EXTRACTION METHOD)  
(GB/T 2912.1-2009 MOD ISO 14104-1:1995)

NOT DETECTED

REQUIREMENT  
≤20 mg/kg  
(Q/AADQ 15-2020)

REMARK: DETECTION LIMIT 20mg/kg

### BACTERIAL FILTRATION EFFICIENCY(%)

(YY 0469-2011 ANNEX B, TEST BACTERIA: STAPHYLOCOCCUS AUREUS ATCC 6508 , TEST AREA:  
40cm<sup>2</sup> , FLOW RATE: 20.0L/min, MEAN PARTICLE SIZE: 3.0μm, RESULT OF THE POSITIVE  
CONTROL:  $2.9 \times 10^3$  CFU , RESULT OF THE NEGATIVE CONTROL: <1CFU )

EFF<sub>1</sub> 99.7  
EFF<sub>2</sub> 99.6  
EFF<sub>3</sub> 99.6

REQUIREMENT  
≥99.5  
(Q/AADQ 15-2020)

### PARTICLE FILTRATION EFFICIENCY(%)

(YY 0469-2011 5.6.2, AIR FLOW:50L/min, AEROSOL:NaCl, AEROSOL CONCENTRATION:15mg/m<sup>3</sup> ,  
TEMP:23.1℃, RH:96.0%)

MINIMUM 97.51

REQUIREMENT  
≥90  
(Q/AADQ 15-2020)

### INSPIRATORY RESISTANCE(Pa)

(GB/T 32610-2016 6.7, HEAD SIZE: MEDIUM)

UNTREATED SAMPLE:

1# 79.9  
2# 81.6

PRETREATMENT SAMPLE:

1# 76.4  
2# 73.9

REQUIREMENT  
≤175  
(Q/AADQ 15-2020)

### EXPIRATORY RESISTANCE(Pa)

(GB/T 32610-2016 6.8, HEAD SIZE: MEDIUM)

UNTREATED SAMPLE:

1# 80.4  
2# 60.6

PRETREATMENT SAMPLE:

1# 80.7  
2# 61.2

REQUIREMENT  
≤185  
(Q/AADQ 15-2020)

### APPEARANCE QUALITY[3 FIBERS]

(Q/AADQ 15-2020 5.1)

PASS

REQUIREMENT要求  
ACCORDING TO THE CLAUSE 4.1 OF THE  
PRODUCT STANDARD  
(Q/AADQ 15-2020)

—End of Report—

PAGE 5 OF 5



# Export approval of disposable mask

## 海关进出口货物收发货人备案回执

企业名称	中天(中国)工业有限公司
统一社会信用代码	91350583777549414P
海关备案日期	2007-01-25
海关编码	3505945316
检验检疫备案号	3501601226
有效期	长期

泉州海关驻南安办事处



2020年06月19日


自然人、法人或者非法人组织可通过“中国海关企业进出口信用信息公示平台”  
(<http://credit.customs.gov.cn>) 或者“互联网+海关”(<http://online.customs.gov.cn>)  
查询海关公示的企业信息。





# AAB口罩

## AAB Mask categories



- 一次性民用防护口罩  
Disposable protective mask

- 一次性使用医用口罩  
Disposable medical mask

- 一次性医用外科口罩  
Disposable surgical mask


- 儿童平面口罩  
Disposable protective mask for children



- KN95防护口罩 (FFP2)  
KN95 protective mask(Equivalent to FFP2)

- KN95防护口罩 (FFP3)  
KN95 protective mask(Equivalent to FFP3)

- KN95医用口罩 (FFP3)  
KN95 medical mask(Equivalent to FFP3)



- 3D 成人防护口罩  
3D Disposable protective mask for adult

- 3D 儿童防护口罩  
3D Disposable protective mask for children



**THANK YOU!**