

Product Specification Sheet

Product description:

- Synthetic Vinyl Powder Free Exam Gloves

Product stock/reorder number:

- 600-V111-CS
- 600-V112-CS
- 600-V113-CS
- 600-V114-CS

Physical characteristics/sizing:

- Small, Medium, Large and X-Large
- Overall length: 240 +/- 5 mm
- Cuff thickness: 0.06 mm
- Palm thickness: 0.08mm
- Finger thickness: 0.05 mm
- Weight: 4.5 grams

Physical properties:

Before aging

- Tensile strength: 15 MPA
- Ultimate elongation: 300%

Composition:

- Polyresin vinyl, DOP (DINP), PVC heat calcium, zinc liquid stabilizers, epoxidized soybean oil

Product features:

- Latex free
- Disposable
- Easy donning
- Powder free
- Ambidextrous

Storage:

- Store in a dry area. Keep away from fire and flame.

Packaging/shipping:

- 10 boxes per case (100 gloves per box)



Professional Medical, Inc.

A Tradition of Quality, Value, & Trusted Service Since 1968

800.648.5190 ■ 1917 Garnet Ct. ■ New Lenox, IL 60451 ■ promedsupply.com