

DIGITAL THERMOMETER (WITH FLEXIBLE TIP)



General Description

A 60-second waterproof digital thermometer with a flexible tip for oral, rectal, or axillary use. Comes with storage case and one battery.

Product Codes

BV Medical
40-130-240-L

Packaging

Unit Packaging (Retail)
Single thermometer with 5 probe covers and case.

Box Packaging
12 ea.

Case Packaging
288 ea.

Specifications

Range: (32.0-42.0)°C / (89.6-107.6)°F
Accuracy: +/-0.1°C/0.2°F (35.5°C-42.0°C/95.9°F-107.6°F), +/-0.2°C/0.4°F (for the remaining measurement range)
Weight: approx. 0.42 oz

BATTERY TYPE

Alkaline manganese battery, type LR41, 1.5V, service life minimum 100 hours under continuous operation.

Features & Benefits

60 SECOND READOUT

A 60 second sensing probe with LCD display with easy to read temperature readings and auto-off feature.

MEMORY RECALL

Able to recall last temperature reading.

PEAK TEMPERATURE

Unit beeps when maximum temperature is reached.

ONE YEAR LIMITED WARRANTY

Barcodes



BV Medical Unit Device Identifier (UDI) / Lot Number

CONTAINS NO MERCURY