

SCHOOL HEALTH DIGITAL THERMOMETER (WATERPROOF)



SPEC SHEET

General Description

A 30 second waterproof digital thermometer with a flexible tip usable for oral, rectal, and axillary. Comes with plastic carrying case and one battery.

Product Code

School Health
#13101

Packaging

Unit Packaging (Retail)
Single thermometer with
5 probe covers.

Box Packaging
12 ea.

Case Packaging
288 ea.

Specifications

Range: (32.0-42.9)°C / (89.6-109.2)°F
Accuracy: +/-0.1°C/0.2°F (35.5°C-42.0°C/95.9°F-107.6°F), +/-0.2°C/0.4°F (32.0°C-35.5°C, 42.0°C-42.9°C/89.6°F-95.9°F, 107.6°F/-109.2°F)
Dimensions: approx. 5 x 0.71 x 0.4 in.
Weight: approx. 0.35 oz

BATTERY TYPE

Alkaline manganese battery, type LR41, 1.5V, service life minimum 100 hours under continuous operation.

Features & Benefits

30 SECOND READOUT

A 30 second sensing probe with LCD display with easy to read temperature readings and auto-off feature.

MEMORY RECALL


Able to recall last temperature reading.

FEVER ALARM

Unit beeps when maximum temperature is reached.

ONE YEAR LIMITED WARRANTY

Barcodes



School Health Reorder #13101



(01)00845821027498(10)0720

Unit Device Identifier (UDI) / Lot Number