

MI 24 & MI 26

Fast, Reliable Tympanometers



Whether you need a reliable tympanometer—the MI 24, or the convenience of a tympanometer with an audiometer—the MI 26, Maico offers this series of tympanometers with the performance, speed and functionality you have come to trust.

- Four frequency, four level reflex capabilities
- Ipsi and contra reflex testing
- Graphical reflex display
- Immediate gradient pass/fail results
- Internal printer and backlit LCD screen
- MI 24 has an optional 1 kHz high frequency probe tone
- MI 26 provides screening air conduction audiometry

EXCELLENCE IN HEARING SCREENING

SINCE 1937



MI 24 & MI 26

Fast, Reliable Tympanometers



Top performance and speed

You demand performance and speed from your tympanometer. Maico provides test repeatability and reliability to meet your most stringent requirements. Results are immediately available on the LCD screen and can be printed automatically on our quiet, fast thermal printer.

Comprehensive testing

The MI 24 and MI 26 perform automatic impedance tests along with multi-frequency, multi-level reflex screening. Both instruments offer graphical reflex display and gradient calculations as well as the option for contralateral reflex testing with either a button phone or TDH 39 headphones. In addition, the MI 24 has a 1 kHz high frequency probe tone option which is ideal for providing reliable results when testing newborns and young children.

Additional testing opportunities

Is office space an issue? With the MI 26 you can bring complete air conduction audiometric testing and tympanometry to your office in one compact instrument. Out-of-office testing is made easier, too, when carrying only one lightweight instrument in our optional soft-side shoulder bag.

CPT codes

MI24 Tympanometer		MI26 Tympanometer	
92567	Tympanometry	92567	Tympanometry
92568	Acoustic reflex	92568	Acoustic reflex
		92551	Screening, pure tone
		92552	Threshold, air

Specifications subject to change without notice

ONE YEAR LIMITED WARRANTY: This warranty is extended to the original purchaser of the instrument, by Maico, through the Distributor from whom it was purchased. The warranty covers defects in material and workmanship for a period of one year from date of delivery of the instrument to the original purchaser. Accessories which are purchased from Maico at the same time as the instrument are warranted for one year from the date of purchase. For additional information contact your Maico Distributor.

EXCELLENCE IN HEARING SCREENING

SINCE 1937



MAICO DIAGNOSTICS • 7625 Golden Triangle Drive, Eden Prairie, MN 55344 • phone 952-941-4200 • toll free 888-941-4201 • fax 952-903-4200 • www.maico-diagnostics.com