

Titmus V3 General Model Vision Screener

Slide Package (SH # 52783)

Slide	Slide Name	Testing Function	Testing Function Detail	Image
1	SAS-PT "E" Symbol Acuity	Acuity, Right & Left Eye	Tests visual acuity (sharpness of the eye). Right & left eye images are different. Administer test at far & using plus lenses.	
2	SPS-3 Muscle Imbalance	Muscle Balance	This is a muscle balance test to test Vertical & Lateral Phoria in children. For children that fail this test, it helps identify if there is a possibility of amblyopia (lazy eye) or strabismus (squint eyes). Administer first at far and then at near.	
3	52-194 LEA SYMBOLS®	Binocularity & Acuity, Left & Right Eye	LEA SYMBOLS® to test Binocularity and Acuity, left & right eye.	
4	SA-BRL-4 Visual Acuity	Acuity, Both, Right & Left Eye	This test checks for visual acuity of both eyes, right eye and left eye using the same slide.	
5	SPCD-1 Pediatric Color	Color Perception	Tests color perception. Tests uses photographically reproduced Titmus Pseudo-Isochromatic blocks.	
6	SPS-4 Combined Phoria	Muscle Balance	This is a Muscle Balance Test to test Vertical and Lateral Phoria in children. For children that FAIL this test, it helps identify if there is a possibility of amblyopia (lazy eye) or strabismus (squint eyes). Test is administered first at far and then at near.	
7	SA-BRL-2A Visual Acuity	Acuity, Both, Right & Left Eye	This test checks for visual acuity of Both Eyes, Right Eye and Left Eye using the same slide	
8	SCD-1 Ishihara Color Perception	Color Perception	Tests red/green color perception. Test uses accurately reproduced Ishihara Pseudo-Isochromatic plates.	