

Titmus V2 General Model Vision Screener

Slide Package (SH # 52789)

| Slide | Slide Name | Testing Function | Testing Function Detail | Image |
|-------|------------------------------------|---|---|-------|
| 1 | SAS-PT "E" Symbol Acuity | Acuity, Right & Left Eye | Tests visual acuity (sharpness of the eye). Right & left eye images are different. Administer test at far & using plus lenses. | |
| 2 | SPS-3 Muscle Imbalance | Muscle Balance | This is a muscle balance test to test Vertical & Lateral Phoria in children. For children that fail this test, it helps identify if there is a possibility of amblyopia (lazy eye) or strabismus (squint eyes). Administer first at far and then at near. | |
| 3 | 52-194 LEA SYMBOLS® | Binocularity & Acuity, Left & Right Eye | LEA SYMBOLS® to test Binocularity and Acuity, left & right eye. | |
| 4 | SA-BRL-4 Visual Acuity | Acuity, Both, Right & Left Eye | This test checks for visual acuity of both eyes, right eye and left eye using the same slide. | |
| 5 | SPCD-1 Pediatric Color | Color Perception | Tests color perception. Tests uses photographically reproduced Titmus Pseudo-Isochromatic blocks. | |
| 6 | SPS-4 Combined Phoria | Muscle Balance | This is a Muscle Balance Test to test Vertical and Lateral Phoria in children. For children that FAIL this test, it helps identify if there is a possibility of amblyopia (lazy eye) or strabismus (squint eyes). Test is administered first at far and then at near. | |
| 7 | SA-BRL-2A Visual Acuity | Acuity, Both, Right & Left Eye | This test checks for visual acuity of Both Eyes, Right Eye and Left Eye using the same slide | |
| 8 | SCD-1 Ishihara Color Perception | Color Perception | Tests red/green color perception. Test uses accurately reproduced Ishihara Pseudo-Isochromatic plates. | |