








Titmus V4 Pediatric Model Vision Screener

Slide Package (SH # 52790)

Slide	Slide Name	Testing Function	Testing Function Detail	Image
1	52-194 LEA SYMBOLS®	Binocularity & Acuity, Left & Right Eye	LEA SYMBOLS® to test Binocularity and Acuity, left & right eye.	 20/100 20/40 20/32 20/25 20/20
2	SMPS-30-1 MI Preschool Acuity	Acuity, Right & Left Eye	Visual acuity testing using the Michigan Preschool Test. Right and Left Eye Images are different. Administer test at FAR.	
3	SMPS-30-2 MI Preschool Acuity	Acuity, Right & Left Eye	Visual acuity testing using the Michigan Preschool Test. Right and Left Eye images are different. Administer test at FAR.	
4	SMPS-30-3 MI Preschool Acuity	Acuity, Right & Left Eye	Visual acuity testing using the Michigan Preschool Test. Right and Left Eye images are different. Administer test at FAR.	
5	SMPS-30-4 MI Preschool Acuity	Acuity, Right & Left Eye	Visual acuity testing using the Michigan Preschool Test. Right and Left Eye images are different. Administer test at FAR.	
6	SAS-PT "E" Symbol Acuity	Acuity, Right & Left Eye	Tests visual acuity (sharpness of the eye). Right and Left Eye images are different. Administer test at FAR and using PLUS lenses. Refer to page 6 for Plus Lens Test.	
7	SPS-3 Muscle Imbalance	Muscle Balance	This is a Muscle Balance Test to test Vertical and Lateral Phoria in children. For children that fail this test, it helps identify if there is a possibility of amblyopia (lazy eye) or strabismus (squint eyes). Test is administered first at far and then at near.	
8	SPCD-1 Pediatric Color	Color Perception	Tests color perception. Tests uses photographically reproduced Titmus Pseudo-Isochromatic blocks.	