

1 Minute Efficacy

Mycobacterium

- *Mycobacterium tuberculosis* var: *bovis* (BCG) (TB)

Bacteria

- *Burkholderia cepacia*
- *Enterobacter cloacae*
- *Klebsiella pneumoniae*
- *Pseudomonas aeruginosa*
- *Salmonella enterica*
- *Staphylococcus aureus*

Pathogenic Fungi

- *Candida albicans*

Drug-resistant Bacteria

- Carbapenem-Resistant *Klebsiella pneumoniae* (CRKP) (CRE)
- ESBL *Escherichia coli* (ESBL E. coli)
- Methicillin Resistant *Staphylococcus aureus* (MRSA)
- Methicillin Resistant *Staphylococcus epidermidis* (MRSE)
- Multi-drug resistant (MDR) *Acinetobacter baumannii*
- Vancomycin Intermediate *Staphylococcus aureus* (VISA)
- Vancomycin Resistant *Enterococcus faecalis* (VRE)

Enveloped Viruses

- Hepatitis B Virus (HBV)
- Hepatitis C Virus (HCV)
- Herpes Simplex Virus Type 1 and Type 2
- Human Coronavirus (not associated with Severe Acute Respiratory Syndrome or SARS)
- Human Immunodeficiency Virus (HIV-1)
- Influenza A Virus Strain A (H3N2 Virus)
- Measles Virus
- SARS-CoV-2 (COVID-19 Virus)

3 Minute Efficacy

Non-Enveloped Viruses

- Adenovirus Type 5*

Fungicide

- *Trichophyton interdigitale**

**Trichophyton interdigitale* requires precleaning. *Trichophyton interdigitale* and Adenovirus require 3-minute contact time.