

## Features:

- Paper dispenser and paper cutter.
- Latex-free.
- 2" thick urethane foam top allows for maximum patient support and comfort. Foam conforms to patient then rebounds to original shape.
- All hardwood leg construction with wood veneer.




Shown in Turquoise with Walnut Base



Shown in Dark Blue with Natural Base



Shown in Black with Natural Base



Shown in Teal with Natural Base

## These Tables Are Built For Athletes.

Heavy-duty construction will make these treatment tables the MVP of your athletic training room.

Feature enclosed storage, keeping supplies dust-free and out of sight.

When ordering, please indicate base color and upholstery color, and laminate color when applicable.

**Handsome:** Your choice of four stain options applied to a beautiful hardwood veneer. Specify Natural or Walnut.

**Tough:** Your choice of 26 durable Texhyde® upholstery colors or 2 high-pressure laminate colors.

**Painless:** Shipped completely assembled on sturdy pallets.

**Practical:** Paper Dispenser and Paper Cutter are standard.

**Safe:** All furnishings are completely latex-free. Security locks come standard.

**Sturdy:** High quality, solid hardwood leg construction (not laminate, composite or particle board).



Royal Blue

Black

Our Exclusive EXPRESS Ship Program will get your SH Treatment Table (D24191) shipped within three days of order receipt. The Three-Day EXPRESS Table is limited to Black or Royal Blue on a Natural base. Other colors (pg. 166) are still available for this treatment table with normal delivery times.

See page 166 for more upholstery and base selections.

## SH Professional Treatment Table with Open Storage

All tables are 78" long, 30" high and have a full-length shelf. Specify your choices from available customization options below.

Appearance	Available Customizations			SH Item #	Price
	Backrest	Face Cutout	Two Fixed Shelves		
A. with 24" Width	✓	✓	✓	*D24212	\$499.00 ea.
with 30" Width	✓	✓	✓	*D24191	\$572.00 ea.
B. with 24" Width		✓	✓	*D24213	\$627.00 ea.
with 30" Width		✓	✓	*D24192	\$699.00 ea.
C. with 30" Width	✓	✓		*D24185	\$616.00 ea.
with 30" Width and Backrest		✓		*D24186	\$754.00 ea.
D. with 24" Width	✓		✓	*D24320	\$578.00 ea.
with 30" Width	✓		✓	*D24321	\$638.00 ea.
with 24" Width and Backrest			✓	*D24322	\$677.00 ea.
with 30" Width and Backrest			✓	*D24324	\$754.00 ea.



## SH Sports Health® Mat Table

Features the same premium hardwood construction as our standard treatment tables. With a 2" foam top upholstered in your choice of Premier Texhyde® vinyl. When ordering, indicate base and upholstery colors. Custom sizes available.

Dimensions: 24"L x 45"W x 21"H

..... \*D24215 \$861.00 ea.

## 6-Legged Treatment Table

This unique design features six legs to distribute weight more easily, and a dropped shelf for extra stability.

Dimensions: 78"L x 30"W x 30"H

with Shelf. .... \*D24497 \$659.00 ea.

with Shelf and Backrest (Shown)

..... \*D24498 \$745.00 ea.

\* Direct-Ship Item. Please specify colors when ordering; allow 3-6 weeks for fabrication and delivery after all necessary information has been submitted. See page 180 for details.



## SH Treatment/Taping Tables & Stations



F



H



E



G

### SH Deluxe Professional Treatment Table with Enclosed Storage

Description	Length	Width	Height	SH Item #	Price
<b>E. with Enclosed Cabinet, Sliding Doors and One Fixed-Position Shelf</b>					
(Shown)	78"	30"	30"	*D24315	\$1,095.00 ea.
Same as D24315 with Backrest	78"	30"	30"	*D24316	\$1,195.00 ea.
<b>F. with Enclosed Central Cabinet and Two Shelves</b>					
(Shown)	78"	30"	30"	*D24187	\$765.00 ea.
Same as D24187 with Backrest	78"	30"	30"	*D24188	\$853.00 ea.
<b>G. with Three 20" Wide Drawers on Nylon Rollers with Positive Stops</b>					
(Shown)	78"	30"	30"	*D24317	\$924.00 ea.
Same as D24317 with Backrest	78"	30"	30"	*D24318	\$999.00 ea.
<b>H. with Three Drawers and Dual Leg Elevation</b>					
(Shown)	78"	30"	30"	*D24332	\$999.00 ea.
Same as D24332 with Backrest	78"	30"	30"	*D24333	\$1,175.00 ea.



#### Full-Length Shelf with Backrest

Features a full-length shelf and adjustable backrest. When ordering, indicate base and upholstery colors. Custom sizes available.

Dimensions: 78" L x 30" H

24" Width	*D24213	\$627.00 ea.
30" Width	*D24192	\$699.00 ea.

\*Direct-Ship Item. Please allow additional time to receive. See page 180 for details.

#### Full-Length Shelf with 2 Fixed Shelves


Features a full-length shelf and two fixed shelves. When ordering, indicate base and upholstery colors. Custom sizes available.

Dimensions: 78" L x 30" H


30" Width	*D24185	\$616.00 ea.
with Adjustable Backrest		
	*D24186	\$754.00 ea.

### How To Order Treatment Tables

#### 1. Select Base Veneer Finish



Natural



Walnut

Please Note: Colors are applied to a wood veneer and will differ slightly in appearance from a laminate of the same color.

#### 2. Select Upholstery

Texhyde®

Leather-Look with Textile Softness



Royal Blue



Parchment



Walnut



Turquoise



Ocean Gray



Med. Green



Dark Green



Teal



Wine



Beige



Light Blue



Bisque



Plum



Dark Gray



Copenhagen



Dark Blue



Black



Stonehedge



Yellow



Jade



Dark Purple



Red



Wedgewood



Apricot



Orange



Tobasco



See page 166 for more upholstery and base selections.

## Dual Taping Station

Dimensions: 69" L x 36" W x 38" H  
Shipping Weight: 460 lbs.

with Upholstered Top ..... \*D24339 \$1,726.00 ea.

with Laminate Top ..... \*D24338 \$1,726.00 ea.

## Single Taping Station (Not Shown)

Dimensions: 44" L x 36" W x 38" H  
Shipping Weight: 318 lbs.

with Upholstered Top ..... \*D24341 \$1,118.00 ea.

with Laminate Top ..... \*D24340 \$1,118.00 ea.



## Taping Table with Two Storage Shelves

Dimensions: 48" L x 24" W x 36" H  
Shipping Weight: 119 lbs.

with Upholstered Top ..... \*D24195 \$463.00 ea.

with Laminate Top ..... \*D24346 \$438.00 ea.

## Taping Table with H-Brace (Not Shown)

Dimensions: 48" L x 24" W x 36" H  
Shipping Weight: 106 lbs.

with Upholstered Top ..... \*D24194 \$356.00 ea.

with Laminate Top ..... \*D24345 \$340.00 ea.



## ProFex Economy Tote-Table

Compact, lightweight and portable. Corner guards and 1" foam top. Folds, making it easy to carry.

Weight: 37 lbs.

Size: 72"L x 24"W x 30"H

Profex Economy Tote-Table, Black ..... \*D24035 \$575.00 ea.

Optional Vinyl Carrying Case, Black ..... \*D24077 \$72.10 ea.



## Portable Massage Table

Made of lightweight aluminum. Easy to store and transport. Covered in oil-resistant, easy to clean, black Boltaflex®.

Standard Height: 30" H

Adjustable Height: 24" H to 30" H

Weight Capacity: 350 lbs.

Standard ..... \*D24347 \$349.00 ea.

Adjustable ..... \*D24348 \$395.00 ea.

Case ..... \*D24350 \$87.00 ea.



## Oakworks Boss Table

Features sealed seams, water-resistant undercoating and puncture-resistant surface which stands up to even the most extreme field conditions. Please specify table width and color when placing order.

Available in 27", 29", and 31" Widths

Weight Capacity: 600 lbs.

Height Range: 24"-34", Length: 73"

Table ..... \*D24352 \$539.00 ea.

Carrying Case ..... \*D24353 \$93.00 ea.

Table & Case ..... \*D24354 \$623.00 ea.

\*Direct-Ship Item. Please allow additional time to receive. See page 180 for details.

## Oakworks Portable Taping Table

Features independently adjustable legs and unique field feet. Marine-grade polymer top and shoulder strap on table.

Color: Black Only

Weight Capacity: 500 lbs., Width: 24"

Height Range: 32"-42", Length: 48"

Table ..... \*D24355 \$494.00 ea.

Carrying Case ..... \*D24356 \$115.00 ea.

Table & Case ..... \*D24357 \$599.00 ea.

## Color Chart for Oakworks Tables



Black



Coffee Cream



Royal Blue



Spruce

\*Direct-Ship Item. Please allow additional time to receive. See page 180 for details.





## Exam Tables

### Ritter® 204 Manual Exam Table

18-gauge steel shell and reinforced seat to support patient, also prevents breakage. Comfort designed stirrups help make the lithotomy position a little easier

er for your patients. Features multi-density foam top design and powder-coated paint finish.

*Dusty Blue, Moss, Navy, Perfect Plum, Shadow*

..... \*D24360 \$1,295.00 ea.



### Patient Positioning System

- Offers superior comfort while positioning patients properly
- Larger vent

areas for easier breathing and advanced lower profile design reduces excessive cervical extension/flexion

- Resilient platform provides optimal support for both prone and supine positioning
- Easily cleanable and durable polymer platform
- Specify Black or Blue when ordering.

..... \*D24359 \$79.00 ea.

### EarthLite® Colors

Upholstery coverings for EarthLite® tables feel soft and warm to the touch. Covering has superior abrasion resistance, is oil- and stain-resistant, and can be cleaned with alcohol or strong detergents.



### EarthLite® Treatment Table

Standard features include EarthLite's special 3½" thick, three-layer foam system, rounded corners for better access and working ease. Headrest outlets are at both ends. Features a hardwood shelf.

*Weight: 78 lbs., Size: 73"L x 30"W*

*Adjustable Height Range: 23"-33"*

*Treatment Table (Specify EarthLite® Color)*

..... \*D24013 \$875.00 ea.

*with Optional Tilt Back (Specify EarthLite® Color)*

..... \*D24014 \$959.00 ea.

\*Direct-Ship Item. Please allow additional time to receive. See page 180 for details.



### Oakworks ProLux PT300

Aluminum base and Tuffles upholstery. Crescent Face Pad can be reversed to support both prone and supine positioning. Please specify table width and color when placing order.

*Black, Coffee Cream, Royal Blue, Spruce*

*Available in 27", 29", and 31" Widths*

*(31" width is an additional \$80.00)*

*Weight Capacity: 550 lbs.*

*Height Range: 18"-34", Length: 73"*

..... \*D24358 \$2,275.00 ea.

\*Direct-Ship Item. Please allow additional time to receive. See page 180 for details.



### EarthLite® Portable Exam Table

Lightweight, yet strong enough to hold more than 2,300 lbs. It's everything you expect from a portable professional table, at an affordable price. High-density, 2" thick foam. Headrests are optional and outlets are at both ends.

*Weight: 30 lbs.*

*Size: 73"L x 30"W*

*Adjustable Height Range: 23"-33"*

*Portable Exam Table with Adjustable Headrest and*

*Black Carrying Case (Specify EarthLite® Color)*

..... \*D24006 \$459.00 ea.

*Adjustable Headrest (Specify EarthLite® Color)*

..... \*D24009 \$79.30 ea.

*Black Carrying Case . . \*D24010 \$106.00 ea.*



### Dual-Wedge Leg and Knee Elevator

High-density urethane foam. Two, vinyl-covered foam wedges with Velcro® brand fasteners. Multifunctional wedge system assembles as a cube or in a variety of configurations.

..... \*D24305 \$139.95 ea.



### Flat Pillow

14" x 14"

.. \*D24054 \$22.00 ea.



### Chatanooga Pron Pillow®

Elevates head and shoulders for added comfort. Raised base allows patients in prone position to breathe easier.

21" x 18" x 6"

..... D24066 \$122.00 ea.

### Paper Roll Dispenser

*with Plastic Band Cutter (for mounting under head)*

..... D24049 \$24.00 ea.

*Plastic Band Cutter Only*

..... D24095 \$5.00 ea.

### Fantastik® All Purpose Cleaner

Perfect for Twill, Texhyde® and laminates.


32 oz. .... D34195 \$4.70 ea.



### Wedge Pillow

25" x 25" x 7", tapered to 1"

.. \*D24106 \$49.00 ea.



### Full- and Half-Round Pillows

Pillows are used as bolsters under the knees, ankles or neck to relieve muscle strain. Made from high-density urethane foam.

*Half-Round Pillow, 3" x 6" x 26"*

..... \*D24291 \$31.00 ea.

*Full-Round Pillow, 6" x 26"*

..... \*D24290 \$31.00 ea.

### Vinyl Liquid Patch

1 oz. .... D90068 \$3.95 ea.

## Medical Cabinetry from SH Sports Health®

Natural wood veneer with the same premium wood construction as our popular SH Sports Health® Treatment Table line. Our cabinetry system is designed to complement other Sports Health® furnishings, making your athletic training room both functional and beautiful. Superior craftsmanship and quality are built into every detail.

\* Direct-Ship Item. Allow additional time to receive. See page 180 for details.

### Standard Color Chart

Please indicate on your order: desired color of cabinet and countertop from these options.

#### Cabinet Color



Natural



Acorn



Auburn



Walnut

#### Countertop Color



White



Black



Fog Dust

#### Drawer/Door Pulls



Metal Wire

All SH Cabinetry features metal pulls.



### Configuration #1

**G** Base Cabinet with **D** Wall Cabinet, shown in Natural.

All Cabinetry features metal pulls and an installed locking system. This system allows the user to unlock all drawers with one key.



### Configuration #2

**H** and **N** Base Cabinets with two **E** Wall Cabinets, shown in Natural.

All multiple configurations consist of heavy-duty interior hardware that attaches units in exact formation. Seamless, one-piece custom countertop with 3 1/2" continuous backsplash.



### Configuration #3

**M**, **K** and **G** Base Cabinets with **B** and **D** Wall Cabinets, shown in Walnut.

When ordering, please specify whether cabinets are to be flush to left wall, right wall or centered with both ends exposed (as shown in these three example configurations).

## Wall Cabinets – 30" H x 12 1/4" D



### 19 1/2"W - Wall Cabinets

**A.** Acrylic..... \*D24233 \$443.00 ea.  
**D.** Solid..... \*D24236 \$346.00 ea.



### 25 1/2"W - Wall Cabinets

**B.** Acrylic..... \*D24234 \$557.00 ea.  
**E.** Solid..... \*D24237 \$428.00 ea.



### 31 1/2"W - Wall Cabinets

**C.** Acrylic..... \*D24235 \$629.00 ea.  
**F.** Solid..... \*D24238 \$485.00 ea.

## Base Cabinets – 36" H x 19" D



### 19 1/2"W - Base Cabinets

**G.** 5 Drawers..... \*D24240 \$768.00 ea.  
**J.** 1 Drawer, 1 Door \*D24245 \$711.00 ea.  
**M.** 1 Door..... \*D24239 \$572.00 ea.



### 25 1/2"W - Base Cabinets

**H.** 5 Drawers..... \*D24243 \$866.00 ea.  
**K.** 2 Drawers, 2 Doors..... \*D24241 \$819.00 ea.  
**N.** 2 Doors..... \*D24246 \$670.00 ea.



### 31 1/2"W - Base Cabinets

**I.** 5 Drawers..... \*D24244 \$979.00 ea.  
**L.** 3 Drawers, 2 Doors..... \*D24242 \$974.00 ea.  
**O.** 2 Doors..... \*D24247 \$716.00 ea.  
**P.** Wheelchair Accessible Sink Cabinet..... \*D24230 \$649.00 ea.



# PROTEAM™

by HAUSMANN

## Proteam™ Exam Tables

### Taping Station Features:

- Upholstered Taping Modules are 36"H or 42"H.
- You can order individual taping or storage modules.
- Durable dowel and glue laminated wood construction.
- Self-closing laminate drawer with nylon rollers.
- Doors have concealed, adjustable, self-closing 110° Euro-hinges.
- All doors and drawer fronts have 3mm high impact edging.
- 2" thick High-Density foam tops are removable without tools.
- Weight capacity of up to 450 lbs.

D



A



B

### Select Series Stations

Description	Length	Height	Depth	SH Item #	Price
<b>A. with Two Seats and Hamper (D24514)</b>	72"	32"	32"	<b>D24499</b>	\$2,275.00
<b>B. with Two Seats</b>	72"	36"	36"	<b>D24500</b>	\$2,175.00
	72"	36"	42"	<b>D24501</b>	\$2,295.00
	72"	36"	48"	<b>D24502</b>	\$2,395.00
	72"	42"	36"	<b>D24503</b>	\$2,295.00
	72"	42"	42"	<b>D24504</b>	\$2,495.00
	72"	42"	48"	<b>D24505</b>	\$2,595.00
<b>C. with Three Seats (Not shown)</b>	120"	36"	36"	<b>D24506</b>	\$3,495.00
<b>D. with Four Seats</b>	168"	36"	36"	<b>D24507</b>	\$4,695.00

Note: D24515 not available on 36" D tables

### Add Custom Logo


Contact Sports Health Customer Service to add a custom logo to your table.



Debossed



Embossed



### Select Custom Options (for Tables B, C, D only)

<b>A. Extra (1) Drawer</b> ...	<b>D24508</b>	\$99.00 ea.
<b>B. Full Extension Pull-Out</b>	<b>D24509</b>	\$86.00 ea.
<b>C. Hinged Side Access Door</b>		
Specify Right or Left...	<b>D24510</b>	\$102.00 ea.
<b>D. 4" Satin Finish Wire Pull</b>		
Free Option	<b>D24511</b>	
<b>E. Velcro Wall Pad, 24" W x 24" D</b>	<b>D24512</b>	\$133.00 ea.

<b>F. 3 Drawers (replacing drawers)</b>	<b>D24513</b>	\$296.00 ea.
<b>G. Waste Hamper and Cut-Out</b>	<b>D24514</b>	\$126.00 ea.
<b>H. Fixed 75° Upholstered Back, 26"H x 10"D</b>	<b>D24515</b>	\$299.99 ea.
<b>I. Laminate Top with 3 mm Edging</b>	<b>D24516</b>	
Free Option		
<b>J. Lock</b>	<b>D24517</b>	\$43.00 ea.
<b>Lock with Bolt</b>	<b>D24518</b>	\$76.00 ea.

### Select Laminate



\* Direct-Ship Item. Please specify color when ordering; allow 3-6 weeks for fabrication and delivery after all information has been submitted. See page 180 for details.

### Select Vinyl Color





## PROTEAM™ by HAUSMANN



### Split Leg Table

- (6) Hinged doors with interior shelf offer concealed supply storage.
- Each padded leg rest can be adjusted and locked into desired angle up to 45° with spring-loaded lock-in lever.
- Comfortable 2" thick High-Density foam upholstered tops.
- Specify Natural Oak, Wild Cherry or Folkstone Gray laminate.

Weight Capacity: 400 lbs.

78" L x 30" W x 31" H

with Open Cabinet... \*D24519 \$1,245.00 ea.

with Hinged Door Shown

..... \*D24520 \$1,595.00 ea.



### 600 lb. Athletic Trainer's Table

- (6) Legs and plywood top for added strength.
- New enhanced leg hardware system.
- Laminate shelf with (2) long stretchers.
- 30" x 78" upholstered top with 2" High-Density foam.
- Legs have protective black vinyl scuff caps.
- Specify Natural Oak or Wild Cherry laminate.

Weight Capacity: 600 lbs.

78" L x 30" W x 31" H

..... \*D24522 \$847.00 ea.



### H-Brace Taping Table

- 36" is convenient height for taping.
- H-Brace is recessed into table leg with unique 4-sided "Lock-Tite" joint.
- High pressure laminate legs, apron and stretchers.
- Legs have protective black vinyl scuff caps.
- 2" thick High-Density foam vinyl upholstered top.
- Only available in Natural Oak laminate.

Weight Capacity: 350 lbs.

48" L x 24" W x 36" H

..... \*D24523 \$385.00 ea.



### Two-Person Taping Table Convertible Treatment Table

- (2) Padded taping sections 26" W x 34" D x 36" H with (1) removable middle Top Pad.
- (3) Storage sections with Oak laminate interiors and shelf.
- (1) Center drawer on nylon rollers.
- Doors have concealed, adjustable, self-closing Euro-hinges.
- 2" High-Density urethane foam vinyl top.
- Only available in Natural Oak laminate.

Weight Capacity: 500 lbs.

74" L x 34" W x 36" H

..... \*D24521 \$1,795.00 ea.



### Pro-Lockers

- Dowel and glue laminated wood construction with 3mm high-impact edging.
- Foot Locker with vented storage and hinged laminate seat (upholstered on D24526) with hasp for padlock.
- Chrome clothes rod on D24524 only.
- Quantity discount available for contract orders.
- Lockers shown with Locking Security Door (for D24524).

32" x 22" x 78" .... \*D24524 \$860.00 ea.

18" x 24" x 78" .... \*D24526 \$110.00 ea.

\*Direct-Ship Item. Please allow additional time to receive. See page 180 for details.