



LIFEPAK CR[®] Plus DEFIBRILLATOR



At Physio-Control, our reputation—and our mission—are built on confidence:

- The confidence that comes from over 50 years of innovation, a steadfast commitment to quality, and a position as the global leader in defibrillation
- The confidence that comes from manufacturing the lifesaving equipment of choice for more EMS and hospital units around the world than any other brand
- The confidence that comes from knowing the LIFEPAK CR Plus defibrillator in your office, the airport or your child's school can simply, safely and effectively help someone save a life

Featuring the same advanced technology trusted by emergency medical professionals—yet simple to use—the LIFEPAK CR Plus automated external defibrillator (AED) is designed specifically for the first person to respond to a victim of sudden cardiac arrest (SCA). Although not everyone can be saved from sudden cardiac arrest, studies show that early defibrillation can dramatically improve survival rates. AED users should be trained in CPR and use of the AED.

Confidence is knowing you will be ready

- Weekly & monthly self tests plus four visible readiness displays let you know if the CR Plus is ready to do its job
- Simple maintenance with synchronized replacement for CHARGE-PAK™ battery charger & QUIK-PAK™ electrodes

Confidence is a simple two-step process

- Smart, safe device designed to shock only if needed
- Designed specifically for minimally-trained users
- Fully automatic model with 2-step operation*
 - Step 1: Turn on device
 - Step 2: Apply electrode pads

*The CR Plus is also available in a semi-automatic model with a shock button.

ONLY EIGHT-YEAR WARRANTY IN THE INDUSTRY



NO NEED TO PUSH A SHOCK BUTTON

Confidence is knowing it's the brand the pros use.

Fast, smooth transition to EMS: QUIK-PAK electrodes easily plug into advanced LIFEPAK® defibrillator/monitor.

More EMS units around the world use LIFEPAK equipment than any other brand, and the electrodes on the *CR Plus* are fully compatible with all other LIFEPAK defibrillators and monitors. That means a faster transfer and more time spent on lifesaving than equipment.



Use Infant/Child Reduced Energy Defibrillation Electrodes on infants and children less than 8 years of age or less than 55 lbs (25

Confidence is maximizing chances for survival.

The only AED on the market today that offers the capability to shock SCA victims from 200 joules to 360 joules. While the initial shock delivered at 200 joules has been shown to be effective at defibrillating the heart of a majority of SCA victims, some patients are more difficult to defibrillate than others. Simply put, the capacity to shock up to 360J gives patients a better opportunity for a better outcome.

- Only AED on the market with energy to escalate from 200 to 300 to 360 joules
- Provides CPR guidance and defibrillation per AHA Guidelines



Included in the price of the *CR Plus*:

- 2 pair QUIK-PAK™ Electrodes
- 1 CHARGE-PAK™ Battery Charger
- Carry Case
- Ambu® Res-Cue Mask® Kit
- 8-year warranty



For more information, please call us!

SCHOOL HEALTH®
We Supply Your Future™

Order Toll Free Phone 866-323-5465 | Order Online www.schoolhealth.com | Order Toll Free Fax 800-235-1305