

Welch Allyn Vital Signs Promotion

UPGRADE AND SAVE

UPGRADE TO A NEW
QUALIFYING DEVICE* TO
RECEIVE A \$150 REBATE

APRIL 1 - SEPTEMBER 30, 2025

***Connex** Spot Monitor, **Connex** Vital Signs Monitor, **Connex** Integrated Wall System, **Spot Vital Signs** 4400 device.

The **Spot Vital Signs** 4200 has been end of manufactured and is scheduled to go End of Service June 30, 2025.

The **Spot Vital Signs** LXi went End of Service June 30, 2023.

Help minimize downtime by upgrading to a new Welch Allyn unit today.

Contact your School Health representative to learn more or redeem this offer at hillrom.com/legacyvitalspromo2025

866-323-5465 | SchoolHealth.com | @SchoolHealth



COMPARISON OF OUR CURRENT SOLUTIONS VS. LEGACY DEVICES

FEATURE	SPOT VITAL SIGNS 4200 DEVICE (END OF MANUFACTURE)	SPOT VITAL SIGNS LXi DEVICE (END OF SERVICE)	SPOT VITAL SIGNS 4400 DEVICE	CONNEX SPOT MONITOR
Workflows	Spot	Spot	Spot	Spot, Office, and Interval Profiles
BP Averaging	No	No	Yes	Yes
Parameters	Step BP, Temp, SPO ₂	SureBP, Temp, SPO ₂	SureBP, Temp, SPO ₂	SureBP, Temp, SPO ₂ , RRp
SPO ₂ Options	Nellcor, Masimo	Nellcor, Masimo	Nonin	Nellcor, Masimo, Nonin
Temperature Options	SureTemp	SureTemp Plus, Braun PRO 4000	SureTemp Plus	SureTemp Plus, Braun PRO 6000
Custom Parameters/Modifiers	None	Basic (height, weight, pain, RR, BMI)	Basic (height, weight, pain, RR, BMI)	Advanced (20 customizable fields)
Early Warning Scores	No	No	No	Yes
Wireless Connectivity	No	Wi-Fi® technology via Accessory	No	Wi-Fi technology, Bluetooth® technology
Wired Connectivity	Serial Cable	USB, Serial Cable	USB	USB or Ethernet
Barcode Scanner	No	Yes	No	Yes
Scale Connectivity	No	Yes via USB Dongle	No	Yes via USB Dongle
Display	5" Monochrome LCD	5.8" Monochrome LCD	7" Color Touchscreen	7" Color Touchscreen
Memory	1 Record	50 records	50 records	400 records
Battery Type	Li Ion	Li Ion	Li Ion	Li Ion
Charge Time	12 Hours	10 Hours	4 Hours	4 Hours
Run Time	8 Hours	8 Hours	8 Hours	8 Hours or 17 Hours with APM
Mounting Options	Basic Stand, Wall Mount, IV Pole Mount	Basic Stand, Wall Mount	Desktop, Basic Stand, 777 Wall Mount	Desktop, Basic Stand, APM Stand, 777 or GCX Wall Mount
School Health Items	—	—	1036782 & 1036783	1034457, 1040300, 1036785, 1036784

Rx only. For safe and proper use of the products mentioned herein, please refer to the appropriate Operator's Manual or Instructions for Use.

Offer available on eligible purchases made April 1 - September 30, 2025. To be eligible, customer must purchase qualifying Welch Allyn products (School Health items 1036782, 1036783, 1034457, 1040300, 1036785 and 1036784) identified in this document and provide proof of Welch Allyn purchase as indicated on the redemption form for this offer found at the link above. Redemption forms and submissions for other promotions/offers cannot be combined. All information must be submitted via the online redemption process by December 15, 2025. All customers must submit a completed W-9. Promotion valid in U.S. for only primary care medical facilities, Long Term Care Facilities, and Ambulatory Surgery Centers. Distributor stocking orders do not qualify. Please allow 12 weeks for receipt of this rebate check based on terms above. Promotion cannot be combined with any other promotion, offer, coupon, enhanced discount or special pricing arrangement; some exclusions may apply.

2505

Baxter, Connex, Spot Vital Signs, SureTemp, and Welch Allyn are trademarks of Baxter International Inc or its subsidiaries.

866-323-5465
SchoolHealth.com
@SchoolHealth

 **SCHOOL HEALTH®**
We Supply Your Future™