

# MobileAid®

Emergency Response



## MobileAid Trauma First Aid Station (31500)

Wheel this 260-piece Trauma First Aid station directly to the medical emergency – large or small.

- Compact station is organized in Quick-Response modules for rapid response and treatment.
- Extensive First Aid supplies meet OSHA requirements and ANSI standards and include a flexible evacuation stretcher.
- High visibility markings, flashing LED wheels, and an 8-foot collapsible flag make the mobile station easy to spot in any situation.
- Unique dual wheel design makes climbing curbs and stairs easy.
- An integrated seat and gear platform allow patient to sit to be treated or can be used to haul additional equipment.
- Super-durable aerospace class frame construction and heavy-duty weather resistant bag safeguard contents.
- Telescoping handle retracts for compact storage in less than 2 cubic feet.

### MobileAid Trauma First Aid Station Contents (260 pieces):

Quick-Response pockets	Quantity	Quick-Response Module B (cont.)	Quantity
First Aid guide book	1	Cotton tip applicator 3"	12
Writing tablet	1	Hydrocortisone 1% (1gram foil pack)	14
Pen	1	Vinyl gloves (pairs)	3
Vinyl gloves (pairs)	2	Emergency water pouches (4 oz)	3
Scissor - utility shears 7-1/4"	1	Hydrogen peroxide solution (4 oz spray)	1
CPR filtershield adult/child	1	Sting & Bite swab	2
Face mask with eyeshield	2	Eye wash 1 oz sterile	2
Protective gown	2		
Flexible emergency evacuation stretcher	1	<b>Quick-Response Module C</b>	
Disposable incident blanket	2	<b>Wound Dressing and Bandaging</b>	
Splint board with pad (medium 12" x 6")	2	Bandage strip 1" x 3"	16
Instruction and inventory sheet	1	Plastic 2" x 3" patch bandage	10
Waterproof document pouch	1	Woven knuckle bandage	10
First Aid flag and pole (collapsible)	1	Gauze pads 3" x 3" 12 ply	8
		Butterfly closure (medium)	16
<b>Quick-Response Module A</b>		Non-adherent pad 2" x 3"	6
<b>Protection, Instruments, Sanitation</b>		Eye pads with adhesive strips	4
First Aid instruction chart	1	Co-Flex bandage 2" x 5 yards	2
Bandage scissors 5-1/2"	1	Gauze clean-wrap bandage 2"	6
Splinter forceps 4-1/2"	1	2" x 6" gauze bandage	2
Lalley forceps 5-1/2"	1	Compress bandage 4" with telfa pad	1
Flashlight (penlight)	1	Bloodstopper bandage	1
AA batteries	4	Triangular bandage 40" x 40" x 56" safety pins	1
Vinyl gloves (pairs)	2	Waterproof adhesive tape 1" x 5 yards	1
Sunblock SPF 30+ (4 oz)	1	Waterproof adhesive tape 1/2" x 5 yards	2
Insect repellent (2 oz)	1		
Emergency blankets	2	<b>Quick-Response Module D</b>	
Hand sanitizer (4 oz)	1	<b>Large Wound Dressing and Bandaging</b>	
Biohazard bag	12	ABD combine pad 5" x 9"	2
Instant towels in a tube	10	Bloodstopper bandage	1
Aspirin (2 5-grain tablets)	25	ABD pads 8" x 10" sterile	4
		Burn dressing gel-soaked 4" x 4"	2
<b>Quick-Response Module B</b>		Sterile gauze roll 4-1/2" x 4.1 yards	2
<b>Wound Cleaning &amp; Treatment</b>		Waterproof adhesive tape 1" x 5 yards	1
First Aid cream (1 gm foil pack)	10		
Burn gel (1/8 oz foil pack)	6	<b>Quick-Response Module E</b>	
PVP iodine swabs (1/2cc each)	10	<b>Hydration &amp; Cold Application</b>	
Wound cleaner towelette	10	Emergency water pouches (4 oz)	3
		Cold pack unit 5" x 6"	4

Equivalent supply substitutions may be made. Product graphics may vary from images shown.

#### Item #31500

Shipping weight & dimensions: 28 lbs; 15.5" x 19.0" x 13.5"

External dimensions: (including wheels) 12.25"W x 18.5"H x 14.0"D

External dimensions: (excluding wheels) 10.0"W x 18.5"H x 14.0"D

Internal dimensions: 8.5"W x 15.5"H x 12.0"D

Access door dimensions: 8.0"W x 13.0"H

Handle extension: 38"

Product weight: 27 lbs

© 2009 LifeSecure. All Rights Reserved. MobileAid is a registered trademark of LifeSecure. MA0909

Be the first responder!™

Respond immediately.

Respond effectively.

Respond confidently.



MobileAid  
Emergency  
Station



MobileAid Trauma  
First Aid Station

For more information about MobileAid Trauma First Aid and Emergency Stations contact your LifeSecure Emergency Solutions dealer.

LifeSecure  
Emergency Solutions